

# 2021 MegaGen **Digital Library** Guide Book





MegaGen Digital Library Guide Book  
© 2021 MegaGen Implant Co., LTD.

**Contact us**

r2b@r2gate.com  
www.r2gate.com

Made by Simon  
(r2.sales3@imegagen.com)



Our Library

exocad

3shape

Dental wings

	ZrGEN (Ti-Base)	TiGEN (Pre-milled)	EZ Post (Stock-Abt)	ZrGEN (Ti-Base)	TiGEN (Pre-milled)	EZ Post (Stock-Abt)	ZrGEN (Ti-Base)	TiGEN (Pre-milled)	EZ Post (Stock-Abt)
<b>AnyRidge</b> Scan abutment Scan Healing abutment	✓	✓	✓	✓	✓	✓	✓	✓	-
<b>BlueDiamond</b> (AnyRidge Octal) Scan abutment Scan Healing abutment	✓	✓	✓	✓	✓	✓	-	-	-
<b>AnyRidge Extra</b> Scan abutment Scan Healing abutment	✓	✓	✓	✓	✓	✓	-	-	-
<b>AnyOne</b> Scan abutment Scan Healing abutment	✓	✓	✓	✓	✓	✓	✓	✓	-
<b>AnyOne OneStage</b> Scan abutment	✓	✓	✓	✓	✓	✓	-	-	-
<b>AnyOne External</b> Scan abutment	✓	✓	-	✓	✓	-	-	-	-
<b>MiNi</b> Scan abutment	✓	✓	✓	✓	✓	✓	✓	✓	-
<b>ST</b> Scan abutment	✓	✓	✓	✓	✓	✓	✓	✓	-
<b>Multi Unit Abutment</b> N type - Scan abutment S type - EZ Post cylinder	✓	-	✓	✓	✓	-	-	-	-
<b>Octa Level</b> Scan abutment	✓	✓	-	✓	✓	-	-	-	-



Scan abutment



Scan Healing abutment



ZrGEN



TiGEN



EZ Post

Library guide

exocad

3shape

KOR

ENG

KOR

ENG

Library Installation

Manual



Scan Healing abutment

Use method manual



EZ Post

Use method manual

-

-




# How to choose library options in CAD software

## exocad




Select an implant system from the list.



When working with ZrGEN,


1. Select the size (Cuff, Post height).
2. Select the cement space.




When working with TiGEN,

1. Select the size (Pre-milled Diameter).


## 3shape



Select an implant Company from the Category




Select the implant system and type (**ZrGEN, TiGEN**) from the system list.



Select an Size (Cuff, Post height Or Pre-milled Diameter) from the Kit list.


# AnyRidge



\* EZ Post Prosthesis by utilizing the stock abutment like a scan abutment.  
\*\* C-type For Sirona CEREC System


# Scan abutment

## Scan Abutment [AR]



Product info			Library info (Choose below CAD library option)		
Diameter	Lenght	Ref.C	exocad library name	3shape library name	Dental wings library name
Ø4.0	9	AANISR4009T	AnyRidge Custom Abutment (9mm)	AnyRidge_ZrGEN_Custom abutment	AnyRidge Internal
	13	AANISR4013T	AnyRidge Custom Abutment (13mm)	AnyRidge_TiGEN_Custom abutment	

## Scan Healing Abutment [AR]



Product info				Library info (Choose below CAD library option)		
Diameter	Height	Abutment Ref.C	Scan Post Ref.C	exocad library name	3shape library name	Dental wings library name
Ø4.0	4	ARISH4004T	SP4007.MTN	AnyRidge SHA Custom Abutment (Healing Height 4mm)	AnyRidge_ZrGEN_SHA_Custom abutment	
	5	ARISH4005T		AnyRidge SHA Custom Abutment (Healing Height 5mm)		
	7	ARISH4007T		AnyRidge SHA Custom Abutment (Healing Height 7mm)		
Ø5.0	4	ARISH5004T	SP5007.MTN	AnyRidge SHA Custom Abutment (Healing Height 4mm)	ZrGEN4015_Healing Height 4mm ZrGEN4015_Healing Height 5mm ZrGEN4015_Healing Height 7mm	
	5	ARISH5005T		AnyRidge SHA Custom Abutment (Healing Height 5mm)		
	7	ARISH5007T		AnyRidge SHA Custom Abutment (Healing Height 7mm)		
Ø6.0	4	ARISH6004T	SP6007.MTN	AnyRidge SHA Custom Abutment (Healing Height 4mm)	AnyRidge_TiGEN_SHA_Custom abutment	
	5	ARISH6005T		AnyRidge SHA Custom Abutment (Healing Height 5mm)		
	7	ARISH6007T		AnyRidge SHA Custom Abutment (Healing Height 7mm)		
Ø7.0	4	ARISH7004T	SP7007.MTN	AnyRidge SHA Custom Abutment (Healing Height 4mm)	TiGEN1020_Healing Height 4mm TiGEN1020_Healing Height 5mm TiGEN1020_Healing Height 7mm	
	5	ARISH7005T		AnyRidge SHA Custom Abutment (Healing Height 5mm)		
	7	ARISH7007T		AnyRidge SHA Custom Abutment (Healing Height 7mm)		

## Scan Abutment C-type [AR] For use CEREC



Product info				
Diameter	Cuff Height	Post height	Post size	Ref. C
Ø3.9	0.5	4.7	Small	ARICSS3405T
	1.0			ARICSS3410T
	2.0			ARICSS3420T
Ø4.3	0.5	4.7	Small	ARICSS3805T
	1.0			ARICSS3810T
	2.0			ARICSS3820T
Ø5.5	0.5	4.7	Large	ARICSL4505T
	1.0			ARICSL4510T
	2.0			ARICSL4520T

\* Compatible with [Xive Library](#) in **CEREC Or in-Lab** CAD Software.

**ZrGEN**  
Ti-Base



Library info		Product info						
ZrGEN Option In exocad library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C		
System	Kit							
	ZrGEN 4015	Ø4.0	0.6	4.5	Hex	AANIPR4015.MTN		
	ZrGEN 4016			6.0		AANIPR4016.MTN		
	ZrGEN 4018			8.0		AANIPR4018.MTN		
	ZrGEN 4025		1.5	4.5		AANIPR4025.MTN		
	ZrGEN 4026			6.0		AANIPR4026.MTN		
	ZrGEN 4028			8.0		AANIPR4028.MTN		
	ZrGEN 4035		3.0	4.5		AANIPR4035.MTN		
	ZrGEN 4036			6.0		AANIPR4036.MTN		
	ZrGEN 4038			8.0		AANIPR4038.MTN		
	ZrGEN 4045		4.0	4.5		AANIPR4045.MTN		
	ZrGEN 4046			6.0		AANIPR4046.MTN		
	ZrGEN 4048			8.0		AANIPR4048.MTN		
	ZrGEN 4015			Ø4.0		0.6	4.5	Non-Hex
	ZrGEN 4016		6.0				AANIPR4016N.MTN	
	ZrGEN 4018	8.0	AANIPR4018N.MTN					
	ZrGEN 4025	1.5	4.5		AANIPR4025N.MTN			
	ZrGEN 4026		6.0		AANIPR4026N.MTN			
	ZrGEN 4028		8.0		AANIPR4028N.MTN			
	ZrGEN 4035	3.0	4.5		AANIPR4035N.MTN			
	ZrGEN 4036		6.0		AANIPR4036N.MTN			
	ZrGEN 4038		8.0		AANIPR4038N.MTN			
	ZrGEN 4045	4.0	4.5		AANIPR4045N.MTN			
	ZrGEN 4046		6.0		AANIPR4046N.MTN			
	ZrGEN 4048		8.0		AANIPR4048N.MTN			
	ZrGEN 4515		Ø4.5		0.6	4.5	Hex	
	ZrGEN 4516	6.0				AANIPR4516.MTN		
	ZrGEN 4518	8.0		AANIPR4518.MTN				
	ZrGEN 4525	1.5		4.5	AANIPR4525.MTN			
	ZrGEN 4526			6.0	AANIPR4526.MTN			
	ZrGEN 4528			8.0	AANIPR4528.MTN			
	ZrGEN 4535	3.0		4.5	AANIPR4535.MTN			
	ZrGEN 4536			6.0	AANIPR4536.MTN			
	ZrGEN 4538			8.0	AANIPR4538.MTN			
	ZrGEN 4545	4.0		4.5	AANIPR4545.MTN			
	ZrGEN 4546			6.0	AANIPR4546.MTN			
	ZrGEN 4548			8.0	AANIPR4548.MTN			
	ZrGEN 4515			Ø4.5	0.6	4.5		Non-Hex
	ZrGEN 4516	6.0				AANIPR4516N.MTN		
	ZrGEN 4518	8.0	AANIPR4518N.MTN					
	ZrGEN 4525	1.5	4.5		AANIPR4525N.MTN			
	ZrGEN 4526		6.0		AANIPR4526N.MTN			
	ZrGEN 4528		8.0		AANIPR4528N.MTN			
	ZrGEN 4535	3.0	4.5		AANIPR4535N.MTN			
	ZrGEN 4536		6.0		AANIPR4536N.MTN			
	ZrGEN 4538		8.0		AANIPR4538N.MTN			
	ZrGEN 4545	4.0	4.5		AANIPR4545N.MTN			
	ZrGEN 4546		6.0		AANIPR4546N.MTN			
	ZrGEN 4548		8.0		AANIPR4548N.MTN			



\* ZrGEN Library provides Cement gap (ZrGEN with Zirconia Abutment space) 0, 10, 20, 30, 40, 50 µ.





**TiGEN**  
Pre-milled






Library info		Product info			
Pre-milled Option In exocad library		Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit				
<p>AnyRidge Custom Abutment(9mm, 13mm)</p>	TiGEN 1020	10	20	Hex	ARTR1020.MTN
	TiGEN 1220	12	20		ARTR1220.MTN
<p>AnyRidge SHA Custom Abutment (Healing Height 4, 5, 7mm)</p>	TiGEN 1020	10	20	Non-Hex	ARTR1020N.MTN
	TiGEN 1220	12	20		ARTR1220N.MTN



**EZ Post**  
Stock abutment



Library info		Product info				
EZ Post Option In exocad library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
	EZ Post 4025	Ø4.0	2.0	5.5	Hex	AANEPH4025L
	EZ Post 4035		3.0			AANEPH4035L
	EZ Post 4045		4.0			AANEPH4045L
	EZ Post 4055		5.0			AANEPH4055L
	EZ Post 4027		2.0	7.0		AANEPH4027L
	EZ Post 4037		3.0			AANEPH4037L
	EZ Post 4047		4.0			AANEPH4047L
	EZ Post 4057		5.0			AANEPH4057L
	EZ Post 4025	Ø4.0	2.0	5.5	Non-Hex	AANEPN4025L
	EZ Post 4035		3.0			AANEPN4035L
	EZ Post 4045		4.0			AANEPN4045L
	EZ Post 4055		5.0			AANEPN4055L
	EZ Post 4027		2.0	7.0		AANEPN4027L
	EZ Post 4037		3.0			AANEPN4037L
	EZ Post 4047		4.0			AANEPN4047L
	EZ Post 4057		5.0			AANEPN4057L
	EZ Post 5025	Ø5.0	2.0	5.5	Hex	AANEPH5025L
	EZ Post 5035		3.0			AANEPH5035L
	EZ Post 5045		4.0			AANEPH5045L
	EZ Post 5055		5.0			AANEPH5055L
	EZ Post 5027		2.0	7.0		AANEPH5027L
	EZ Post 5037		3.0			AANEPH5037L
	EZ Post 5047		4.0			AANEPH5047L
	EZ Post 5057		5.0			AANEPH5057L
	EZ Post 5025	Ø5.0	2.0	5.5	Non-Hex	AANEPN5025L
	EZ Post 5035		3.0			AANEPN5035L
	EZ Post 5045		4.0			AANEPN5045L
	EZ Post 5055		5.0			AANEPN5055L
	EZ Post 5027		2.0	7.0		AANEPN5027L
	EZ Post 5037		3.0			AANEPN5037L
	EZ Post 5047		4.0			AANEPN5047L
	EZ Post 5057		5.0			AANEPN5057L
	EZ Post 6025	Ø6.0	2.0	5.5	Hex	AANEPH6025L
	EZ Post 6035		3.0			AANEPH6035L
	EZ Post 6045		4.0			AANEPH6045L
	EZ Post 6055		5.0			AANEPH6055L
	EZ Post 6027		2.0	7.0		AANEPH6027L
	EZ Post 6037		3.0			AANEPH6037L
	EZ Post 6047		4.0			AANEPH6047L
	EZ Post 6057		5.0			AANEPH6057L
	EZ Post 6025	Ø6.0	2.0	5.5	Non-Hex	AANEPN6025L
	EZ Post 6035		3.0			AANEPN6035L
	EZ Post 6045		4.0			AANEPN6045L
	EZ Post 6055		5.0			AANEPN6055L
	EZ Post 6027		2.0	7.0		AANEPN6027L
	EZ Post 6037		3.0			AANEPN6037L
	EZ Post 6047		4.0			AANEPN6047L
	EZ Post 6057		5.0			AANEPN6057L
	EZ Post 7025	Ø7.0	2.0	5.5	Hex	AANEPH7025L
	EZ Post 7035		3.0			AANEPH7035L
	EZ Post 7045		4.0			AANEPH7045L
	EZ Post 7055		5.0			AANEPH7055L
	EZ Post 7027		2.0	7.0		AANEPH7027L
	EZ Post 7037		3.0			AANEPH7037L
	EZ Post 7047		4.0			AANEPH7047L
	EZ Post 7057		5.0			AANEPH7057L
	EZ Post 7025	Ø7.0	2.0	5.5	Non-Hex	AANEPN7025L
	EZ Post 7035		3.0			AANEPN7035L
	EZ Post 7045		4.0			AANEPN7045L
	EZ Post 7055		5.0			AANEPN7055L
	EZ Post 7027		2.0	7.0		AANEPN7027L
	EZ Post 7037		3.0			AANEPN7037L
	EZ Post 7047		4.0			AANEPN7047L
	EZ Post 7057		5.0			AANEPN7057L



**ZrGEN**  
Ti-Base



Library info			Product info						
ZrGEN Option In 3shape library			Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C		
System	Kit for 9mm	Kit for 13mm							
	ZrGEN 4015_S09	ZrGEN 4015_S13	Ø4.0	0.6	4.5	Hex	AANIPR4015.MTN		
	ZrGEN 4016_S09	ZrGEN 4016_S13			6.0		AANIPR4016.MTN		
	ZrGEN 4018_S09	ZrGEN 4018_S13			8.0		AANIPR4018.MTN		
	ZrGEN 4025_S09	ZrGEN 4025_S13		1.5	4.5		AANIPR4025.MTN		
	ZrGEN 4026_S09	ZrGEN 4026_S13			6.0		AANIPR4026.MTN		
	ZrGEN 4028_S09	ZrGEN 4028_S13			8.0		AANIPR4028.MTN		
	ZrGEN 4035_S09	ZrGEN 4035_S13		3.0	4.5		AANIPR4035.MTN		
	ZrGEN 4036_S09	ZrGEN 4036_S13			6.0		AANIPR4036.MTN		
	ZrGEN 4038_S09	ZrGEN 4038_S13			8.0		AANIPR4038.MTN		
	ZrGEN 4045_S09	ZrGEN 4045_S13		4.0	4.5		AANIPR4045.MTN		
	ZrGEN 4046_S09	ZrGEN 4046_S13			6.0		AANIPR4046.MTN		
	ZrGEN 4048_S09	ZrGEN 4048_S13			8.0		AANIPR4048.MTN		
	ZrGEN 4015_S09	ZrGEN 4015_S13		Ø4.0	0.6		4.5	Non-Hex	AANIPR4015N.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13					6.0		AANIPR4016N.MTN
	ZrGEN 4018_S09	ZrGEN 4018_S13	8.0			AANIPR4018N.MTN			
	ZrGEN 4025_S09	ZrGEN 4025_S13	1.5		4.5	AANIPR4025N.MTN			
	ZrGEN 4026_S09	ZrGEN 4026_S13			6.0	AANIPR4026N.MTN			
	ZrGEN 4028_S09	ZrGEN 4028_S13			8.0	AANIPR4028N.MTN			
	ZrGEN 4035_S09	ZrGEN 4035_S13	3.0		4.5	AANIPR4035N.MTN			
	ZrGEN 4036_S09	ZrGEN 4036_S13			6.0	AANIPR4036N.MTN			
	ZrGEN 4038_S09	ZrGEN 4038_S13			8.0	AANIPR4038N.MTN			
	ZrGEN 4045_S09	ZrGEN 4045_S13	4.0		4.5	AANIPR4045N.MTN			
	ZrGEN 4046_S09	ZrGEN 4046_S13			6.0	AANIPR4046N.MTN			
	ZrGEN 4048_S09	ZrGEN 4048_S13			8.0	AANIPR4048N.MTN			
	ZrGEN 4515_S09	ZrGEN 4515_S13	Ø4.5		0.6	4.5	Hex		AANIPR4515.MTN
	ZrGEN 4516_S09	ZrGEN 4516_S13				6.0			AANIPR4516.MTN
	ZrGEN 4518_S09	ZrGEN 4518_S13		8.0		AANIPR4518.MTN			
	ZrGEN 4525_S09	ZrGEN 4525_S13		1.5	4.5	AANIPR4525.MTN			
	ZrGEN 4526_S09	ZrGEN 4526_S13			6.0	AANIPR4526.MTN			
	ZrGEN 4528_S09	ZrGEN 4528_S13			8.0	AANIPR4528.MTN			
	ZrGEN 4535_S09	ZrGEN 4535_S13		3.0	4.5	AANIPR4535.MTN			
	ZrGEN 4536_S09	ZrGEN 4536_S13			6.0	AANIPR4536.MTN			
	ZrGEN 4538_S09	ZrGEN 4538_S13			8.0	AANIPR4538.MTN			
	ZrGEN 4545_S09	ZrGEN 4545_S13		4.0	4.5	AANIPR4545.MTN			
	ZrGEN 4546_S09	ZrGEN 4546_S13			6.0	AANIPR4546.MTN			
	ZrGEN 4548_S09	ZrGEN 4548_S13			8.0	AANIPR4548.MTN			
	ZrGEN 4515_S09	ZrGEN 4515_S13		Ø4.5	0.6	4.5		Non-Hex	AANIPR4515N.MTN
	ZrGEN 4516_S09	ZrGEN 4516_S13				6.0			AANIPR4516N.MTN
	ZrGEN 4518_S09	ZrGEN 4518_S13	8.0			AANIPR4518N.MTN			
	ZrGEN 4525_S09	ZrGEN 4525_S13	1.5		4.5	AANIPR4525N.MTN			
	ZrGEN 4526_S09	ZrGEN 4526_S13			6.0	AANIPR4526N.MTN			
	ZrGEN 4528_S09	ZrGEN 4528_S13			8.0	AANIPR4528N.MTN			
	ZrGEN 4535_S09	ZrGEN 4535_S13	3.0		4.5	AANIPR4535N.MTN			
	ZrGEN 4536_S09	ZrGEN 4536_S13			6.0	AANIPR4536N.MTN			
	ZrGEN 4538_S09	ZrGEN 4538_S13			8.0	AANIPR4538N.MTN			
	ZrGEN 4545_S09	ZrGEN 4545_S13	4.0		4.5	AANIPR4545N.MTN			
	ZrGEN 4546_S09	ZrGEN 4546_S13			6.0	AANIPR4546N.MTN			
	ZrGEN 4548_S09	ZrGEN 4548_S13			8.0	AANIPR4548N.MTN			



  
AnyRidge\_ZrGEN\_  
Custom abutment

  
AnyRidge\_ZrGEN\_SHA\_  
Custom abutment



**TiGEN**  
Pre-milled






Library info			Product info			
Pre-milled Option In 3shape library			Diameter (∅)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm				
 AnyRidge_TiGEN_ Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	10	20	Hex	ARTR1020.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	12	20		ARTR1220.MTN
 AnyRidge_TiGEN_SHA Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	10	20	Non-Hex	ARTR1020N.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	12	20		ARTR1220N.MTN



**EZ Post**  
Stock abutment



Library info		Product info				
EZ Post Option In 3shape library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
	EZ Post 4025	Ø4.0	2.0	5.5	Hex	AANEPH4025L
	EZ Post 4035		3.0			AANEPH4035L
	EZ Post 4045		4.0			AANEPH4045L
	EZ Post 4055		5.0			AANEPH4055L
	EZ Post 4027		2.0	7.0		AANEPH4027L
	EZ Post 4037		3.0			AANEPH4037L
	EZ Post 4047		4.0			AANEPH4047L
	EZ Post 4057		5.0			AANEPH4057L
	EZ Post 4025	Ø4.0	2.0	5.5	Non-Hex	AANEPN4025L
	EZ Post 4035		3.0			AANEPN4035L
	EZ Post 4045		4.0			AANEPN4045L
	EZ Post 4055		5.0			AANEPN4055L
	EZ Post 4027		2.0	7.0		AANEPN4027L
	EZ Post 4037		3.0			AANEPN4037L
	EZ Post 4047		4.0			AANEPN4047L
	EZ Post 4057		5.0			AANEPN4057L
 AnyRidge_EZ Post_Scan abutment	EZ Post 5025	Ø5.0	2.0	5.5	Hex	AANEPH5025L
	EZ Post 5035		3.0			AANEPH5035L
	EZ Post 5045		4.0			AANEPH5045L
	EZ Post 5055		5.0			AANEPH5055L
	EZ Post 5027		2.0	7.0		AANEPH5027L
	EZ Post 5037		3.0			AANEPH5037L
	EZ Post 5047		4.0			AANEPH5047L
	EZ Post 5057		5.0			AANEPH5057L
 AnyRidge_EZ Post_SHA	EZ Post 5025	Ø5.0	2.0	5.5	Non-Hex	AANEPN5025L
	EZ Post 5035		3.0			AANEPN5035L
	EZ Post 5045		4.0			AANEPN5045L
	EZ Post 5055		5.0			AANEPN5055L
	EZ Post 5027		2.0	7.0		AANEPN5027L
	EZ Post 5037		3.0			AANEPN5037L
	EZ Post 5047		4.0			AANEPN5047L
	EZ Post 5057		5.0			AANEPN5057L
 AnyRidge_EZ Post_EZ Post	EZ Post 6025	Ø6.0	2.0	5.5	Hex	AANEPH6025L
	EZ Post 6035		3.0			AANEPH6035L
	EZ Post 6045		4.0			AANEPH6045L
	EZ Post 6055		5.0			AANEPH6055L
	EZ Post 6027		2.0	7.0		AANEPH6027L
	EZ Post 6037		3.0			AANEPH6037L
	EZ Post 6047		4.0			AANEPH6047L
	EZ Post 6057		5.0			AANEPH6057L
	EZ Post 6025	Ø6.0	2.0	5.5	Non-Hex	AANEPN6025L
	EZ Post 6035		3.0			AANEPN6035L
	EZ Post 6045		4.0			AANEPN6045L
	EZ Post 6055		5.0			AANEPN6055L
	EZ Post 6027		2.0	7.0		AANEPN6027L
	EZ Post 6037		3.0			AANEPN6037L
	EZ Post 6047		4.0			AANEPN6047L
	EZ Post 6057		5.0			AANEPN6057L
	EZ Post 7025	Ø7.0	2.0	5.5	Hex	AANEPH7025L
	EZ Post 7035		3.0			AANEPH7035L
	EZ Post 7045		4.0			AANEPH7045L
	EZ Post 7055		5.0			AANEPH7055L
	EZ Post 7027		2.0	7.0		AANEPH7027L
	EZ Post 7037		3.0			AANEPH7037L
	EZ Post 7047		4.0			AANEPH7047L
	EZ Post 7057		5.0			AANEPH7057L
	EZ Post 7025	Ø7.0	2.0	5.5	Non-Hex	AANEPN7025L
	EZ Post 7035		3.0			AANEPN7035L
	EZ Post 7045		4.0			AANEPN7045L
	EZ Post 7055		5.0			AANEPN7055L
	EZ Post 7027		2.0	7.0		AANEPN7027L
	EZ Post 7037		3.0			AANEPN7037L
	EZ Post 7047		4.0			AANEPN7047L
	EZ Post 7057		5.0			AANEPN7057L

**ZrGEN**

Ti-Base




Library info			Product info				
ZrGEN Option In Dental wings library			Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
	ZrGEN 4015_S09	ZrGEN 4015_S13	Ø4.0	0.6	4.5	Hex	AANIPR4015.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13			6.0		AANIPR4016.MTN
	ZrGEN 4015_S09	ZrGEN 4015_S13	Ø4.0	0.6	4.5	Non-Hex	AANIPR4015N.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13			6.0		AANIPR4016N.MTN
AnyRidge Internal ZrGEN (A_07A6)	ZrGEN 4525_S09	ZrGEN 4525_S13	Ø4.5	0.6	4.5	Hex	AANIPR4525.MTN
	ZrGEN 4526_S09	ZrGEN 4526_S13			6.0		AANIPR45266.MTN
	ZrGEN 4525_S09	ZrGEN 4525_S13	Ø4.5	0.6	4.5	Non-Hex	AANIPR4525N.MTN
	ZrGEN 4526_S09	ZrGEN 4526_S13			6.0		AANIPR4526N.MTN



AnyRidge Internal  
ZrGEN (A\_07A6)





Library info			Product info			
Pre-milled Option In Dental wings library			Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm				
 AnyRidge Internal ✓ TiGEN (A_07A6)	TiGEN 1220_S09	TiGEN 1220_S13	12	20	Hex	ARTR1220.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	12	20	Non-Hex	ARTR1220N.MTN



**ZrGEN**  
For CEREC




Product info				
Diameter (Ø)	Cuff height (mm)	Post height (mm)	Post size (mm)	Ref. C
Ø3.9	0.5	4.7	Small	ARICSS3405.MTN
	1.0			ARICSS3410.MTN
	2.0			ARICSS3420.MTN
Ø4.3	0.5			ARICSS3805.MTN
	1.0			ARICSS3810.MTN
	2.0			ARICSS3820.MTN
Ø5.5	0.5		Large	ARICSL4505.MTN
	1.0			ARICSL4510.MTN
	2.0			ARICSL4520.MTN

\* Compatible with [Xive Library](#) in **CEREC Or in-Lab** CAD Software.





# BlueDiamond



\* EZ Post Prosthesis by utilizing the stock abutment like a scan abutment.  
 \*\* C-type For Sirona CEREC System


# Scan abutment

## Scan Abutment [BD]



Product info				Library info (Choose below CAD library option)		
Diameter	Lenght	Type	Ref.C	exocad library name	3shape library name	Dental wings library name
Ø4.0	13	NC	AROSANT	AnyRidge Octa1_Narrow Custom Abutment (13mm)	AnyRidge Octa1_Narrow_ZrGEN_Custom abutment	-
	13	RC	AROSART	AnyRidge Octa1_Regular Custom Abutment (13mm)	AnyRidge Octa1_Narrow_TiGEN_Custom abutment	-

## Scan Healing Abutment [BD]



Product info					Library info (Choose below CAD library option)		
Type	Diameter (D)	Height	Abutment Ref.C	Scan Post Ref.C	exocad library name	3shape library name	Dental wings library name
NC	Ø4.0	4	AROISHN4004T	SP4007.MTN	AnyRidge Octa1 SHA Narrow Custom Abutment (Healing Height 4mm)	For ZrGEN AnyRidge Octa1_Narrow_ZrGEN_SHA_Custom abutment AnyRidge Octa1_Regular_ZrGEN_SHA_Custom abutment	-
		5	AROISHN4005T		AnyRidge Octa1 SHA Narrow Custom Abutment (Healing Height 5mm)		
		7	AROISHN4007T		AnyRidge Octa1 SHA Narrow Custom Abutment (Healing Height 7mm)		
	Ø5.0	4	AROISHN5004T	SP5007.MTN	AnyRidge Octa1 SHA Narrow Custom Abutment (Healing Height 4mm)		
		5	AROISHN5005T		AnyRidge Octa1 SHA Narrow Custom Abutment (Healing Height 5mm)		
		7	AROISHN5007T		AnyRidge Octa1 SHA Narrow Custom Abutment (Healing Height 7mm)		
RC	Ø4.0	4	AROISHR4004T	SP4007.MTN	AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 4mm)	ZrGEN4015_Healing Height 4mm ZrGEN4015_Healing Height 5mm ZrGEN4015_Healing Height 7mm	-
		5	AROISHR4005T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 5mm)		
		7	AROISHR4007T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 7mm)		
	Ø5.0	4	AROISHR5007T	SP5007.MTN	AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 4mm)	For TiGEN	
		5	AROISHR5005T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 5mm)		
		7	AROISHR5007T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 7mm)		
	Ø6.0	4	AROISHR6004T	SP6007.MTN	AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 4mm)	AnyRidge Octa1_TiGEN_SHA_Custom abutment	
		5	AROISHR6005T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 5mm)		
		7	AROISHR6007T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 7mm)		
	Ø7.0	4	AROISHR7004T	SP7007.MTN	AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 4mm)	TiGEN1020_Healing Height 4mm TiGEN1020_Healing Height 5mm TiGEN1020_Healing Height 7mm	
		5	AROISHR7005T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 5mm)		
		7	AROISHR7007T		AnyRidge Octa1 SHA Regular Custom Abutment (Healing Height 7mm)		

## Scan Abutment C-type [BD] For use CEREC



Product info							
Type	Diameter (D)	Cuff Height	Post height	Post size	Ref. C		
NC	Ø3.9	0.5	4.7	Small	AROCSS3405NT		
		1.0			AROCSS3410NT		
		2.0			AROCSS3420NT		
	Ø4.3	0.5			AROCSS3805NT		
		1.0			AROCSS3810NT		
		2.0			AROCSS3820NT		
RC	Ø3.9	0.5	4.7	Small	AROCSS3405RT		
		1.0			AROCSS3410RT		
		2.0			AROCSS3420RT		
	Ø4.3	0.5			AROCSS3805RT		
		1.0			AROCSS3810RT		
		2.0			AROCSS3820RT		
	Ø5.5	0.5			Large	AROCSS4505RT	
		1.0				AROCSS4510RT	
		2.0				AROCSS4510RT	

\* Compatible with [Xive Library](#) in **CEREC Or in-Lab** CAD Software.

**ZrGEN**  
Ti-Base



AnyRidge Octal  
Narrow Custom Abutment







AnyRidge Octal SHA  
Narrow Custom Abutment  
(Healing Height 4, 5, 7mm)



AnyRidge Octal  
Regular Custom Abutment



AnyRidge Octal SHA  
Regular Custom Abutment  
(Healing Height 4, 5, 7mm)

Library info		Product info					
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit						
 AnyRidge Octal Narrow Custom Abutment   AnyRidge Octal SHA Narrow Custom Abutment (Healing Height 4, 5, 7mm)	ZrGEN 4015	NC	Ø4.0	0.6	4.5	Octa	AROZGN4015.MTN
	ZrGEN 4016				6.0		AROZGN4025.MTN
	ZrGEN 4018				8.0		AROZGN4035.MTN
	ZrGEN 4025			1.5	4.5		AROZGN4045.MTN
	ZrGEN 4026				6.0		AROZGN4016.MTN
	ZrGEN 4028				8.0		AROZGN4026.MTN
	ZrGEN 4035			3.0	4.5		AROZGN4036.MTN
	ZrGEN 4036				6.0		AROZGN4046.MTN
	ZrGEN 4038				8.0		AROZGN4018.MTN
	ZrGEN 4045			4.0	4.5		AROZGN4028.MTN
	ZrGEN 4046				6.0		AROZGN4038.MTN
	ZrGEN 4048				8.0		AROZGN4048.MTN
 AnyRidge Octal Regular Custom Abutment   AnyRidge Octal SHA Regular Custom Abutment (Healing Height 4, 5, 7mm)	ZrGEN 4015	RC	Ø4.5	0.6	4.5	Octa	AROZGN4515.MTN
	ZrGEN 4016				6.0		AROZGN4525.MTN
	ZrGEN 4018				8.0		AROZGN4535.MTN
	ZrGEN 4025			1.5	4.5		AROZGN4535.MTN
	ZrGEN 4026				6.0		AROZGN4516.MTN
	ZrGEN 4028				8.0		AROZGN4526.MTN
	ZrGEN 4035			3.0	4.5		AROZGN4536.MTN
	ZrGEN 4036				6.0		AROZGN4546.MTN
	ZrGEN 4038				8.0		AROZGN4518.MTN
	ZrGEN 4045			4.0	4.5		AROZGN4528.MTN
	ZrGEN 4046				6.0		AROZGN4538.MTN
	ZrGEN 4048				8.0		AROZGN4548.MTN

\* ZrGEN Library provides Cement gap (ZrGEN and Abutment space) 0, 10, 20, 30, 40, 50 µ.





Pre-milled














Library info		Product info				
Pre-milled Option In exocad library		Size	Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit					
 AnyRidge Octa1_TiGEN _Custom abutment	TiGEN 1020	NC	10	20	Octa	AROTGN1020.MTN
	TiGEN 1220		12	20		AROTGN1220.MTN
 AnyRidge Octa1_TiGEN_SHA _Custom abutment	TiGEN 1020	RC	10	20	Octa	AROTGR1020.MTN
	TiGEN 1220		12	20		AROTGR1220.MTN



**EZ Post**  
Stock abutment



Library info		Product info													
EZ Post Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C								
System	Kit														
 AnyRidge Octal_Narrow Scanbody - EZpost (13mm)	EZ Post 4014	NC	Ø4.0	1.0	4.0	Octa	AROEPN4014T								
	EZ Post 4024			2.0			AROEPN4024T								
	EZ Post 4034			3.0			AROEPN4034T								
	EZ Post 4044			4.0			AROEPN4044T								
	EZ Post 4054			5.0			AROEPN4054T								
	EZ Post 4015			1.0	5.5		AROEPN4015T								
	EZ Post 4025			2.0			AROEPN4025T								
	EZ Post 4035			3.0			AROEPN4035T								
	EZ Post 4045			4.0			AROEPN4045T								
	EZ Post 4055			5.0			AROEPN4055T								
	 AnyRidge Octal_Narrow SHA - EZpost (Healing Height 4, 5, 7mm)			EZ Post 4017	NC		Ø4.0	1.0	7.0	Octa	AROEPN4017T				
				EZ Post 4027				2.0			AROEPN4027T				
				EZ Post 4037				3.0			AROEPN4037T				
				EZ Post 4047				4.0			AROEPN4047T				
				EZ Post 4057				5.0			AROEPN4057T				
			 AnyRidge Octal_Narrow EZpost - EZpost	EZ Post 5014		NC		Ø5.0	1.0		4.0	Octa	AROEPN5014T		
				EZ Post 5024					2.0				AROEPN5024T		
				EZ Post 5034					3.0				AROEPN5034T		
				EZ Post 5044					4.0				AROEPN5044T		
				EZ Post 5054					5.0				AROEPN5054T		
	 AnyRidge Octal_Regular Scanbody - EZpost (13mm)			EZ Post 5015	RC		Ø5.0		1.0	5.5	Octa		AROEPN5015T		
				EZ Post 5025					2.0				AROEPN5025T		
				EZ Post 5035					3.0				AROEPN5035T		
				EZ Post 5045					4.0				AROEPN5045T		
				EZ Post 5055					5.0				AROEPN5055T		
 AnyRidge Octal_Regular SHA - EZpost (Healing Height 4, 5, 7mm)		EZ Post 5017	RC	Ø5.0		1.0		7.0	Octa	AROEPN5017T					
		EZ Post 5027				2.0				AROEPN5027T					
		EZ Post 5037				3.0				AROEPN5037T					
		EZ Post 5047				4.0				AROEPN5047T					
		EZ Post 5057				5.0				AROEPN5057T					
		 AnyRidge Octal_Regular EZpost - EZpost				EZ Post 6014		RC		Ø6.0		1.0	4.0	Octa	AROEPN6014T
						EZ Post 6024						2.0			AROEPN6024T
						EZ Post 6034						3.0			AROEPN6034T
						EZ Post 6044						4.0			AROEPN6044T
						EZ Post 6054						5.0			AROEPN6054T
 AnyRidge Octal_Regular EZpost - EZpost			EZ Post 6015	RC		Ø6.0	1.0		5.5		Octa	AROEPN6015T			
			EZ Post 6025				2.0					AROEPN6025T			
			EZ Post 6035				3.0					AROEPN6035T			
			EZ Post 6045				4.0					AROEPN6045T			
			EZ Post 6055				5.0					AROEPN6055T			
		 AnyRidge Octal_Regular EZpost - EZpost	EZ Post 6017				RC	Ø6.0	1.0	7.0		Octa	AROEPN6017T		
			EZ Post 6027						2.0				AROEPN6027T		
			EZ Post 6037						3.0				AROEPN6037T		
			EZ Post 6047						4.0				AROEPN6047T		
			EZ Post 6057						5.0				AROEPN6057T		
	 AnyRidge Octal_Regular EZpost - EZpost		EZ Post 7014		RC				Ø7.0	1.0			4.0	Octa	AROEPN7014T
			EZ Post 7024							2.0					AROEPN7024T
			EZ Post 7034							3.0					AROEPN7034T
			EZ Post 7044							4.0					AROEPN7044T
			EZ Post 7054							5.0					AROEPN7054T
 AnyRidge Octal_Regular EZpost - EZpost		EZ Post 7015	RC	Ø7.0		1.0	5.5	Octa		AROEPN7015T					
		EZ Post 7025				2.0				AROEPN7025T					
		EZ Post 7035				3.0				AROEPN7035T					
		EZ Post 7045				4.0				AROEPN7045T					
		EZ Post 7055				5.0				AROEPN7055T					
		 AnyRidge Octal_Regular EZpost - EZpost				EZ Post 7017	RC			Ø7.0	1.0	7.0	Octa		AROEPN7017T
						EZ Post 7027					2.0				AROEPN7027T
						EZ Post 7037					3.0				AROEPN7037T
						EZ Post 7047					4.0				AROEPN7047T
						EZ Post 7057					5.0				AROEPN7057T

**ZrGEN**

Ti-Base



AnyRidge Octal\_Narrow\_ZrGEN\_Custom abutment



AnyRidge Octal\_Narrow\_ZrGEN\_SHA\_Custom abutment



AnyRidge Octal\_Regular\_ZrGEN\_Custom abutment





AnyRidge Octal\_Regular\_ZrGEN\_SHA\_Custom abutment

Library info			Product info					
ZrGEN Option In 3shape library			Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm						
<p>AnyRidge Octal_Narrow_ZrGEN_Custom abutment</p>	ZrGEN 4015_S09	ZrGEN 4015_S13	NC	Ø4.0	0.6	4.5	Octa	AROZGN4015.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13				6.0		AROZGN4025.MTN
	ZrGEN 4018_S09	ZrGEN 4018_S13				8.0		AROZGN4035.MTN
	ZrGEN 4025_S09	ZrGEN 4025_S13			1.5	4.5		AROZGN4045.MTN
	ZrGEN 4026_S09	ZrGEN 4026_S13				6.0		AROZGN4016.MTN
	ZrGEN 4028_S09	ZrGEN 4028_S13				8.0		AROZGN4026.MTN
	ZrGEN 4035_S09	ZrGEN 4035_S13			3.0	4.5		AROZGN4036.MTN
	ZrGEN 4036_S09	ZrGEN 4036_S13				6.0		AROZGN4046.MTN
	ZrGEN 4038_S09	ZrGEN 4038_S13				8.0		AROZGN4018.MTN
	ZrGEN 4045_S09	ZrGEN 4045_S13			4.0	4.5		AROZGN4028.MTN
	ZrGEN 4046_S09	ZrGEN 4046_S13				6.0		AROZGN4038.MTN
	ZrGEN 4048_S09	ZrGEN 4048_S13				8.0		AROZGN4048.MTN
<p>AnyRidge Octal_Regular_ZrGEN_Custom abutment</p>	ZrGEN 4015_S09	ZrGEN 4015_S13	RC	Ø4.5	0.6	4.5	Octa	AROZGN4515.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13				6.0		AROZGN4525.MTN
	ZrGEN 4018_S09	ZrGEN 4018_S13				8.0		AROZGN4535.MTN
	ZrGEN 4025_S09	ZrGEN 4025_S13			1.5	4.5		AROZGN4535.MTN
	ZrGEN 4026_S09	ZrGEN 4026_S13				6.0		AROZGN4516.MTN
	ZrGEN 4028_S09	ZrGEN 4028_S13				8.0		AROZGN4526.MTN
	ZrGEN 4035_S09	ZrGEN 4035_S13			3.0	4.5		AROZGN4536.MTN
	ZrGEN 4036_S09	ZrGEN 4036_S13				6.0		AROZGN4546.MTN
	ZrGEN 4038_S09	ZrGEN 4038_S13				8.0		AROZGN4518.MTN
	ZrGEN 4045_S09	ZrGEN 4045_S13			4.0	4.5		AROZGN4528.MTN
	ZrGEN 4046_S09	ZrGEN 4046_S13				6.0		AROZGN4538.MTN
	ZrGEN 4048_S09	ZrGEN 4048_S13				8.0		AROZGN4548.MTN



**TiGEN**  
Pre-milled










Library info			Product info				
Pre-milled Option In 3shape library			Size	Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
 AnyRidge_Octa1_TiGEN_ Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	NC	10	20	Octa	AROTGN1020.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		AROTGN1220.MTN
 AnyRidge Octa1_TiGEN_ SHA_Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	RC	10	20	Octa	AROTGR1020N.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		AROTGR1220N.MTN




**EZ Post**  
Stock abutment



Library info		Product info										
EZ Post Option In 3shape library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C					
System	Kit											
 <p>AnyRidge Octa1_Narrow_ EZ Post_Scan abutment</p>	EZ Post 4014	NC	Ø4.0	1.0	4.0	Octa	AROEPN4014T					
	EZ Post 4024			2.0			AROEPN4024T					
	EZ Post 4034			3.0			AROEPN4034T					
	EZ Post 4044			4.0			AROEPN4044T					
	EZ Post 4054			5.0			AROEPN4054T					
	EZ Post 4015			1.0	5.5		AROEPN4015T					
	EZ Post 4025			2.0			AROEPN4025T					
	EZ Post 4035			3.0			AROEPN4035T					
	EZ Post 4045			4.0			AROEPN4045T					
	EZ Post 4055			5.0			AROEPN4055T					
	 <p>AnyRidge Octa1_Narrow_ EZ Post_Scan Healing Abutment</p>		EZ Post 4017	NC	Ø4.0		1.0	7.0	AROEPN4017T			
			EZ Post 4027				2.0		AROEPN4027T			
			EZ Post 4037				3.0		AROEPN4037T			
			EZ Post 4047				4.0		AROEPN4047T			
			EZ Post 4057				5.0		AROEPN4057T			
			 <p>AnyRidge Octa1_Narrow_ EZ Post_EZ Post</p>				EZ Post 5014	NC	Ø5.0	1.0	4.0	AROEPN5014T
							EZ Post 5024			2.0		AROEPN5024T
							EZ Post 5034			3.0		AROEPN5034T
							EZ Post 5044			4.0		AROEPN5044T
							EZ Post 5054			5.0		AROEPN5054T
 <p>AnyRidge Octa1_Regular_ EZ Post_Scan abutment</p>	EZ Post 5015	RC		Ø5.0	1.0	5.5	AROEPN5015T					
	EZ Post 5025				2.0		AROEPN5025T					
	EZ Post 5035				3.0		AROEPN5035T					
	EZ Post 5045				4.0		AROEPN5045T					
	EZ Post 5055				5.0		AROEPN5055T					
	 <p>AnyRidge Octa1_Regular_ EZ Post_Scan Healing Abutment</p>		EZ Post 5017		RC	Ø5.0	1.0	7.0	AROEPN5017T			
			EZ Post 5027				2.0		AROEPN5027T			
			EZ Post 5037				3.0		AROEPN5037T			
			EZ Post 5047				4.0		AROEPN5047T			
			EZ Post 5057				5.0		AROEPN5057T			
			 <p>AnyRidge Octa1_Regular_ EZ Post_EZ Post</p>	EZ Post 6014			RC	Ø6.0	1.0	4.0	AROEPN6014T	
				EZ Post 6024					2.0		AROEPN6024T	
				EZ Post 6034					3.0		AROEPN6034T	
				EZ Post 6044					4.0		AROEPN6044T	
				EZ Post 6054					5.0		AROEPN6054T	
	 <p>AnyRidge Octa1_Regular_ EZ Post_EZ Post</p>			EZ Post 6015	RC	Ø6.0		1.0	5.5	AROEPN6015T		
				EZ Post 6025				2.0		AROEPN6025T		
				EZ Post 6035				3.0		AROEPN6035T		
				EZ Post 6045				4.0		AROEPN6045T		
				EZ Post 6055				5.0		AROEPN6055T		
<p>AnyRidge Octa1_Regular_ EZ Post_EZ Post</p>		EZ Post 6017	RC	Ø7.0		1.0	4.0	AROEPN6017T				
		EZ Post 6027				2.0		AROEPN6027T				
		EZ Post 6037				3.0		AROEPN6037T				
		EZ Post 6047				4.0		AROEPN6047T				
		EZ Post 6057				5.0		AROEPN6057T				
	<p>AnyRidge Octa1_Regular_ EZ Post_EZ Post</p>	EZ Post 7014		RC	Ø7.0	1.0	5.5	AROEPN7014T				
		EZ Post 7024				2.0		AROEPN7024T				
		EZ Post 7034				3.0		AROEPN7034T				
		EZ Post 7044				4.0		AROEPN7044T				
		EZ Post 7054				5.0		AROEPN7054T				
<p>AnyRidge Octa1_Regular_ EZ Post_EZ Post</p>		EZ Post 7015	RC		Ø7.0	1.0	7.0	AROEPN7015T				
		EZ Post 7025				2.0		AROEPN7025T				
		EZ Post 7035				3.0		AROEPN7035T				
		EZ Post 7045				4.0		AROEPN7045T				
		EZ Post 7055				5.0		AROEPN7055T				
	<p>AnyRidge Octa1_Regular_ EZ Post_EZ Post</p>	EZ Post 7017		RC	Ø7.0	1.0	4.0	AROEPN7017T				
		EZ Post 7027				2.0		AROEPN7027T				
		EZ Post 7037				3.0		AROEPN7037T				
		EZ Post 7047				4.0		AROEPN7047T				
		EZ Post 7057				5.0		AROEPN7057T				




# AnyRidge Extra



# Scan abutment

**Scan Abutment [AR]** \* Compatible with **AnyRidge** Scan abutment



Product info			Library info (Choose below CAD library option)		
Diameter	Lenght	Ref.C	exocad library name	3shape library name	Dental wings library name
Ø4.0	9	AANISR4009T	AnyRidge Extra Custom Abutment (9mm)	AnyRidge Extra_Core3.3_ZrGEN_Custom abutment Core3.8, Core4.0, Core4.8	-
	13	AANISR4013T	AnyRidge Extra Custom Abutment (13mm)	AnyRidge Extra_TiGEN_Custom abutment	-



**ZrGEN**  
Ti-Base



Library info		Product info					
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit						
Core 3.3	ZrGEN4515	Core 3.3	Ø4.5	0.6	4.5	Hex	ARZXM503815.MTN
	ZrGEN4516				6.0		ARZXM503816.MTN
	ZrGEN4518				8.0		ARZXM503818.MTN
	ZrGEN4525			1.5	4.5		ARZXM503825.MTN
	ZrGEN4526				6.0		ARZXM503826.MTN
	ZrGEN4528				8.0		ARZXM503828.MTN
	ZrGEN4535			3.0	4.5		ARZXM503835.MTN
	ZrGEN4536				6.0		ARZXM503836.MTN
	ZrGEN4538				8.0		ARZXM503838.MTN
	ZrGEN4545			4.0	4.5		ARZXM503845.MTN
	ZrGEN4546				6.0		ARZXM503846.MTN
	ZrGEN4548				8.0		ARZXM503848.MTN
	ZrGEN4515		Ø4.5	0.6	4.5	Non-hex	ARZXM553815.MTN
	ZrGEN4516				6.0		ARZXM553816.MTN
	ZrGEN4518				8.0		ARZXM553818.MTN
	ZrGEN4525			1.5	4.5		ARZXM553825.MTN
	ZrGEN4526				6.0		ARZXM553826.MTN
	ZrGEN4528				8.0		ARZXM553828.MTN
	ZrGEN4535			3.0	4.5		ARZXM553835.MTN
	ZrGEN4536				6.0		ARZXM553836.MTN
	ZrGEN4538				8.0		ARZXM553838.MTN
	ZrGEN4545			4.0	4.5		ARZXM553845.MTN
	ZrGEN4546				6.0		ARZXM553846.MTN
	ZrGEN4548				8.0		ARZXM553848.MTN
Core 3.8	ZrGEN503815	Core 3.8	Ø5.0	0.6	4.5	Hex	ARZXM503815N.MTN
	ZrGEN503816				6.0		ARZXM503816N.MTN
	ZrGEN503818				8.0		ARZXM503818N.MTN
	ZrGEN503825			1.5	4.5		ARZXM503825N.MTN
	ZrGEN503826				6.0		ARZXM503826N.MTN
	ZrGEN503828				8.0		ARZXM503828N.MTN
	ZrGEN503835			3.0	4.5		ARZXM503835N.MTN
	ZrGEN503836				6.0		ARZXM503836N.MTN
	ZrGEN503838				8.0		ARZXM503838N.MTN
	ZrGEN503845			4.0	4.5		ARZXM503845N.MTN
	ZrGEN503846				6.0		ARZXM503846N.MTN
	ZrGEN503848				8.0		ARZXM503848N.MTN
	ZrGEN503815		Ø5.0	0.6	4.5	Non-hex	ARZXM553815N.MTN
	ZrGEN503816				6.0		ARZXM553816N.MTN
	ZrGEN503818				8.0		ARZXM553818N.MTN
	ZrGEN503825			1.5	4.5		ARZXM553825N.MTN
	ZrGEN503826				6.0		ARZXM553826N.MTN
	ZrGEN503828				8.0		ARZXM553828N.MTN
	ZrGEN503835			3.0	4.5		ARZXM553835N.MTN
	ZrGEN503836				6.0		ARZXM553836N.MTN
	ZrGEN503838				8.0		ARZXM553838N.MTN
	ZrGEN503845			4.0	4.5		ARZXM553845N.MTN
	ZrGEN503846				6.0		ARZXM553846N.MTN
	ZrGEN503848				8.0		ARZXM553848N.MTN
Core 3.8	ZrGEN553815	Core 3.8	Ø5.5	0.6	4.5	Hex	ARZXM553815.MTN
	ZrGEN553816				6.0		ARZXM553816.MTN
	ZrGEN553818				8.0		ARZXM553818.MTN
	ZrGEN553825			1.5	4.5		ARZXM553825.MTN
	ZrGEN553826				6.0		ARZXM553826.MTN
	ZrGEN553828				8.0		ARZXM553828.MTN
	ZrGEN553835			3.0	4.5		ARZXM553835.MTN
	ZrGEN553836				6.0		ARZXM553836.MTN
	ZrGEN553838				8.0		ARZXM553838.MTN
	ZrGEN553845			4.0	4.5		ARZXM553845.MTN
	ZrGEN553846				6.0		ARZXM553846.MTN
	ZrGEN553848				8.0		ARZXM553848.MTN
	ZrGEN553815		Ø5.5	0.6	4.5	Non-hex	ARZXM553815N.MTN
	ZrGEN553816				6.0		ARZXM553816N.MTN
	ZrGEN553818				8.0		ARZXM553818N.MTN
	ZrGEN553825			1.5	4.5		ARZXM553825N.MTN
	ZrGEN553826				6.0		ARZXM553826N.MTN
	ZrGEN553828				8.0		ARZXM553828N.MTN
	ZrGEN553835			3.0	4.5		ARZXM553835N.MTN
	ZrGEN553836				6.0		ARZXM553836N.MTN
	ZrGEN553838				8.0		ARZXM553838N.MTN
	ZrGEN553845			4.0	4.5		ARZXM553845N.MTN
	ZrGEN553846				6.0		ARZXM553846N.MTN
	ZrGEN553848				8.0		ARZXM553848N.MTN



**ZrGEN**

Ti-Base



Library info		Product info					
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit						
AnyRidge Extra Custom Abutment	ZrGEN5015	Core 4.0	Ø5.0	0.6	4.5	Hex	ARZXM5015.MTN
	ZrGEN5016				6.0		ARZXM5016.MTN
	ZrGEN5018				8.0		ARZXM5018.MTN
	ZrGEN5025			1.5	4.5		ARZXM5025.MTN
	ZrGEN5026				6.0		ARZXM5026.MTN
	ZrGEN5028				8.0		ARZXM5028.MTN
	ZrGEN5035			3.0	4.5		ARZXM5035.MTN
	ZrGEN5036				6.0		ARZXM5036.MTN
	ZrGEN5038				8.0		ARZXM5038.MTN
	ZrGEN5045			4.0	4.5		ARZXM5045.MTN
	ZrGEN5046				6.0		ARZXM5046.MTN
	ZrGEN5048				8.0		ARZXM5048.MTN
	ZrGEN5015		Ø5.0	0.6	4.5	Non-hex	ARZXM5015N.MTN
	ZrGEN5016				6.0		ARZXM5016N.MTN
	ZrGEN5018				8.0		ARZXM5018N.MTN
	ZrGEN5025			1.5	4.5		ARZXM5025N.MTN
	ZrGEN5026				6.0		ARZXM5026N.MTN
	ZrGEN5028				8.0		ARZXM5028N.MTN
	ZrGEN5035			3.0	4.5		ARZXM5035N.MTN
	ZrGEN5036				6.0		ARZXM5036N.MTN
ZrGEN5038	8.0	ARZXM5038N.MTN					
ZrGEN5045	4.0	4.5		ARZXM5045N.MTN			
ZrGEN5046		6.0		ARZXM5046N.MTN			
ZrGEN5048		8.0		ARZXM5048N.MTN			
ZrGEN5515	Core 4.0	Ø5.5	0.6	4.5	Hex	ARZXM5515.MTN	
ZrGEN5516				6.0		ARZXM5516.MTN	
ZrGEN5518				8.0		ARZXM5518.MTN	
ZrGEN5525			1.5	4.5		ARZXM5525.MTN	
ZrGEN5526				6.0		ARZXM5526.MTN	
ZrGEN5528				8.0		ARZXM5528.MTN	
ZrGEN5535			3.0	4.5		ARZXM5535.MTN	
ZrGEN5536				6.0		ARZXM5536.MTN	
ZrGEN5538				8.0		ARZXM5538.MTN	
ZrGEN5545			4.0	4.5		ARZXM5545.MTN	
ZrGEN5546				6.0		ARZXM5546.MTN	
ZrGEN5548				8.0		ARZXM5548.MTN	
ZrGEN5515		Ø5.5	0.6	4.5	Non-hex	ARZXM5515N.MTN	
ZrGEN5516				6.0		ARZXM5516N.MTN	
ZrGEN5518				8.0		ARZXM5518N.MTN	
ZrGEN5525			1.5	4.5		ARZXM5525N.MTN	
ZrGEN5526				6.0		ARZXM5526N.MTN	
ZrGEN5528				8.0		ARZXM5528N.MTN	
ZrGEN5535			3.0	4.5		ARZXM5535N.MTN	
ZrGEN5536				6.0		ARZXM5536N.MTN	
ZrGEN5538	8.0			ARZXM5538N.MTN			
ZrGEN5545	4.0		4.5	ARZXM5545N.MTN			
ZrGEN5546			6.0	ARZXM5546N.MTN			
ZrGEN5548			8.0	ARZXM5548N.MTN			



AnyRidge Extra Custom Abutment



**ZrGEN**  
Ti-Base



Library info		Product info					
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit						
AnyRidge Extra Custom Abutment	ZrGEN5515	Core 4.8	Ø5.0	0.6	4.5	Hex	ARZXL5515.MTN
	ZrGEN5516				6.0		ARZXL5516.MTN
	ZrGEN5518				8.0		ARZXL5518.MTN
	ZrGEN5525			1.5	4.5		ARZXL5525.MTN
	ZrGEN5526				6.0		ARZXL5526.MTN
	ZrGEN5528				8.0		ARZXL5528.MTN
	ZrGEN5535			3.0	4.5		ARZXL5535.MTN
	ZrGEN5536				6.0		ARZXL5536.MTN
	ZrGEN5538				8.0		ARZXL5538.MTN
	ZrGEN5545			4.0	4.5		ARZXL5545.MTN
	ZrGEN5546				6.0		ARZXL5546.MTN
	ZrGEN5548				8.0		ARZXL5548.MTN
	ZrGEN5515		Ø5.0	0.6	4.5	Non-hex	ARZXL5515N.MTN
	ZrGEN5516				6.0		ARZXL5516N.MTN
	ZrGEN5518				8.0		ARZXL5518N.MTN
	ZrGEN5525			1.5	4.5		ARZXL5525N.MTN
	ZrGEN5526				6.0		ARZXL5526N.MTN
	ZrGEN5528				8.0		ARZXL5528N.MTN
	ZrGEN5535			3.0	4.5		ARZXL5535N.MTN
	ZrGEN5536				6.0		ARZXL5536N.MTN
	ZrGEN5538				8.0		ARZXL5538N.MTN
	ZrGEN5545			4.0	4.5		ARZXL5545N.MTN
	ZrGEN5546				6.0		ARZXL5546N.MTN
	ZrGEN5548				8.0		ARZXL5548N.MTN
ZrGEN6015	Core 4.8	Ø5.5	0.6	4.5	Hex	ARZXL6015.MTN	
ZrGEN6016				6.0		ARZXL6016.MTN	
ZrGEN6018				8.0		ARZXL6018.MTN	
ZrGEN6025			1.5	4.5		ARZXL6025.MTN	
ZrGEN6026				6.0		ARZXL6026.MTN	
ZrGEN6028				8.0		ARZXL6028.MTN	
ZrGEN6035			3.0	4.5		ARZXL6035.MTN	
ZrGEN6036				6.0		ARZXL6036.MTN	
ZrGEN6038				8.0		ARZXL6038.MTN	
ZrGEN6045			4.0	4.5		ARZXL6045.MTN	
ZrGEN6046				6.0		ARZXL6046.MTN	
ZrGEN6048				8.0		ARZXL6048.MTN	
ZrGEN6015		Ø5.5	0.6	4.5	Non-hex	ARZXL6015N.MTN	
ZrGEN6016				6.0		ARZXL6016N.MTN	
ZrGEN6018				8.0		ARZXL6018N.MTN	
ZrGEN6025			1.5	4.5		ARZXL6025N.MTN	
ZrGEN6026				6.0		ARZXL6026N.MTN	
ZrGEN6028				8.0		ARZXL6028N.MTN	
ZrGEN6035			3.0	4.5		ARZXL6035N.MTN	
ZrGEN6036				6.0		ARZXL6036N.MTN	
ZrGEN6038				8.0		ARZXL6038N.MTN	
ZrGEN6045			4.0	4.5		ARZXL6045N.MTN	
ZrGEN6046				6.0		ARZXL6046N.MTN	
ZrGEN6048				8.0		ARZXL6048N.MTN	



Pre-milled



Bevel Platform



AnyRidge Extra Custom Abutment (9, 13mm)

Library info		Product info				
Pre-milled Option In exocad library		Core (∅)	Diameter (∅)	Length (mm)	Type	Ref.C
System	Kit					
AnyRidge Extra Custom Abutment (9, 13mm)	TiGEN 1020	3.3 Bevel	10	20	Hex	ARTXN1020.MTN
	TiGEN 1220		12	20		ARTXN1020N.MTN
	TiGEN 1020		10	20	Non-Hex	ARTXN1220.MTN
	TiGEN 1220		12	20		ARTXN1220N.MTN
	TiGEN 1020	4.0 Platform	10	20	Hex	ARTXM1020.MTN
	TiGEN 1220		12	20		ARTXM1020N.MTN
	TiGEN 1020		10	20	Non-Hex	ARTXM1220.MTN
	TiGEN 1220		12	20		ARTXM1220N.MTN
	TiGEN 1020	4.8 Platform	10	20	Hex	ARTXL1020.MTN
	TiGEN 1220		12	20		ARTXL1020N.MTN
	TiGEN 1020		10	20	Non-Hex	ARTXL1220.MTN
	TiGEN 1220		12	20		ARTXL1220N.MTN



**ZrGEN**

Ti-Base



Library info			Product info								
ZrGEN Option In 3shape library			Core (Ø)	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C			
System	Kit for 9mm	Kit for 13mm									
AnyRidge_Extra_ZrGEN_Custom abutment	ZrGEN4515_S09	ZrGEN4515_S13	Core 3.3	Ø4.5	0.6	4.5	Hex	ARZNX4515.MTN			
	ZrGEN4516_S09	ZrGEN4516_S13				6.0		ARZNX4516.MTN			
	ZrGEN4518_S09	ZrGEN4518_S13				8.0		ARZNX4518.MTN			
	ZrGEN4525_S09	ZrGEN4525_S13			1.5	4.5		ARZNX4525.MTN			
	ZrGEN4526_S09	ZrGEN4526_S13				6.0		ARZNX4526.MTN			
	ZrGEN4528_S09	ZrGEN4528_S13				8.0		ARZNX4528.MTN			
	ZrGEN4535_S09	ZrGEN4535_S13			3.0	4.5		ARZNX4535.MTN			
	ZrGEN4536_S09	ZrGEN4536_S13				6.0		ARZNX4536.MTN			
	ZrGEN4538_S09	ZrGEN4538_S13				8.0		ARZNX4538.MTN			
	ZrGEN4545_S09	ZrGEN4545_S13			4.0	4.5	ARZNX4545.MTN				
	ZrGEN4546_S09	ZrGEN4546_S13				6.0	ARZNX4546.MTN				
	ZrGEN4548_S09	ZrGEN4548_S13				8.0	ARZNX4548.MTN				
	ZrGEN4515_S09	ZrGEN4515_S13			Core 3.3	Ø4.5	0.6	4.5	Non-hex	ARZNX4515.MTN	
	ZrGEN4516_S09	ZrGEN4516_S13						6.0		ARZNX4516.MTN	
	ZrGEN4518_S09	ZrGEN4518_S13						8.0		ARZNX4518.MTN	
	ZrGEN4525_S09	ZrGEN4525_S13					1.5	4.5		ARZNX4525.MTN	
	ZrGEN4526_S09	ZrGEN4526_S13						6.0		ARZNX4526.MTN	
	ZrGEN4528_S09	ZrGEN4528_S13						8.0		ARZNX4528.MTN	
	ZrGEN4535_S09	ZrGEN4535_S13		3.0			4.5	ARZNX4535.MTN			
	ZrGEN4536_S09	ZrGEN4536_S13					6.0	ARZNX4536.MTN			
	ZrGEN4538_S09	ZrGEN4538_S13					8.0	ARZNX4538.MTN			
	ZrGEN4545_S09	ZrGEN4545_S13		4.0			4.5	ARZNX4545.MTN			
	ZrGEN4546_S09	ZrGEN4546_S13					6.0	ARZNX4546.MTN			
	ZrGEN4548_S09	ZrGEN4548_S13					8.0	ARZNX4548.MTN			
	ZrGEN503815_S09	ZrGEN503815_S13		Core 3.8			Ø5.0	0.6	4.5	Hex	ARZXM503815.MTN
	ZrGEN503816_S09	ZrGEN503816_S13							6.0		ARZXM503816.MTN
	ZrGEN503818_S09	ZrGEN503818_S13							8.0		ARZXM503818.MTN
	ZrGEN503825_S09	ZrGEN503825_S13						1.5	4.5		ARZXM503825.MTN
	ZrGEN503826_S09	ZrGEN503826_S13							6.0		ARZXM503826.MTN
	ZrGEN503828_S09	ZrGEN503828_S13							8.0		ARZXM503828.MTN
	ZrGEN503835_S09	ZrGEN503835_S13			3.0	4.5		ARZXM503835.MTN			
	ZrGEN503836_S09	ZrGEN503836_S13				6.0		ARZXM503836.MTN			
	ZrGEN503838_S09	ZrGEN503838_S13				8.0		ARZXM503838.MTN			
	ZrGEN503845_S09	ZrGEN503845_S13			4.0	4.5		ARZXM503845.MTN			
	ZrGEN503846_S09	ZrGEN503846_S13				6.0		ARZXM503846.MTN			
	ZrGEN503848_S09	ZrGEN503848_S13				8.0		ARZXM503848.MTN			
ZrGEN503815_S09	ZrGEN503815_S13	Core 3.8	Ø5.0		0.6	4.5		Non-hex	ARZXM503815.MTN		
ZrGEN503816_S09	ZrGEN503816_S13					6.0			ARZXM503816.MTN		
ZrGEN503818_S09	ZrGEN503818_S13					8.0			ARZXM503818.MTN		
ZrGEN503825_S09	ZrGEN503825_S13				1.5	4.5			ARZXM503825.MTN		
ZrGEN503826_S09	ZrGEN503826_S13					6.0			ARZXM503826.MTN		
ZrGEN503828_S09	ZrGEN503828_S13					8.0			ARZXM503828.MTN		
ZrGEN503835_S09	ZrGEN503835_S13				3.0	4.5	ARZXM503835.MTN				
ZrGEN503836_S09	ZrGEN503836_S13					6.0	ARZXM503836.MTN				
ZrGEN503838_S09	ZrGEN503838_S13					8.0	ARZXM503838.MTN				
ZrGEN503845_S09	ZrGEN503845_S13				4.0	4.5	ARZXM503845.MTN				
ZrGEN503846_S09	ZrGEN503846_S13					6.0	ARZXM503846.MTN				
ZrGEN503848_S09	ZrGEN503848_S13					8.0	ARZXM503848.MTN				
ZrGEN553815_S09	ZrGEN553815_S13				Core 3.8	Ø5.5	0.6	4.5	Hex	ARZXM553815.MTN	
ZrGEN553816_S09	ZrGEN553816_S13							6.0		ARZXM553816.MTN	
ZrGEN553818_S09	ZrGEN553818_S13							8.0		ARZXM553818.MTN	
ZrGEN553825_S09	ZrGEN553825_S13						1.5	4.5		ARZXM553825.MTN	
ZrGEN553826_S09	ZrGEN553826_S13							6.0		ARZXM553826.MTN	
ZrGEN553828_S09	ZrGEN553828_S13							8.0		ARZXM553828.MTN	
ZrGEN553835_S09	ZrGEN553835_S13	3.0	4.5				ARZXM553835.MTN				
ZrGEN553836_S09	ZrGEN553836_S13		6.0				ARZXM553836.MTN				
ZrGEN553838_S09	ZrGEN553838_S13		8.0				ARZXM553838.MTN				
ZrGEN553845_S09	ZrGEN553845_S13	4.0	4.5				ARZXM553845.MTN				
ZrGEN553846_S09	ZrGEN553846_S13		6.0				ARZXM553846.MTN				
ZrGEN553848_S09	ZrGEN553848_S13		8.0				ARZXM553848.MTN				
ZrGEN553815_S09	ZrGEN553815_S13	Core 3.8	Ø5.5	0.6			4.5	Non-hex	ARZXM553815.MTN		
ZrGEN553816_S09	ZrGEN553816_S13						6.0		ARZXM553816.MTN		
ZrGEN553818_S09	ZrGEN553818_S13						8.0		ARZXM553818.MTN		
ZrGEN553825_S09	ZrGEN553825_S13			1.5			4.5		ARZXM553825.MTN		
ZrGEN553826_S09	ZrGEN553826_S13						6.0		ARZXM553826.MTN		
ZrGEN553828_S09	ZrGEN553828_S13						8.0		ARZXM553828.MTN		
ZrGEN553835_S09	ZrGEN553835_S13			3.0		4.5	ARZXM553835.MTN				
ZrGEN553836_S09	ZrGEN553836_S13					6.0	ARZXM553836.MTN				
ZrGEN553838_S09	ZrGEN553838_S13					8.0	ARZXM553838.MTN				
ZrGEN553845_S09	ZrGEN553845_S13			4.0		4.5	ARZXM553845.MTN				
ZrGEN553846_S09	ZrGEN553846_S13					6.0	ARZXM553846.MTN				
ZrGEN553848_S09	ZrGEN553848_S13					8.0	ARZXM553848.MTN				



**ZrGEN**

Ti-Base



Library info			Product info						
ZrGEN Option In 3shape library			Core (Ø)	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C	
System	Kit for 9mm	Kit for 13mm							
	ZrGEN5015_S09	ZrGEN5015_S13	Core 4.8	Ø5.0	0.6	4.5	Hex	ARZXM5015.MTN	
	ZrGEN5016_S09	ZrGEN5016_S13				6		ARZXM5016.MTN	
	ZrGEN5018_S09	ZrGEN5018_S13				8		ARZXM5018.MTN	
	ZrGEN5025_S09	ZrGEN5025_S13			1.5	4.5		ARZXM5025.MTN	
	ZrGEN5026_S09	ZrGEN5026_S13				6		ARZXM5026.MTN	
	ZrGEN5028_S09	ZrGEN5028_S13				8		ARZXM5028.MTN	
	ZrGEN5035_S09	ZrGEN5035_S13			3.0	4.5		ARZXM5035.MTN	
	ZrGEN5036_S09	ZrGEN5036_S13				6		ARZXM5036.MTN	
	ZrGEN5038_S09	ZrGEN5038_S13				8		ARZXM5038.MTN	
	ZrGEN5045_S09	ZrGEN5045_S13			4.0	4.5		ARZXM5045.MTN	
	ZrGEN5046_S09	ZrGEN5046_S13				6		ARZXM5046.MTN	
	ZrGEN5048_S09	ZrGEN5048_S13				8		ARZXM5048.MTN	
	ZrGEN5015_S09	ZrGEN5015_S13		Core 4.8	Ø5.0	0.6	4.5	Non-hex	ARZXM5015N.MTN
	ZrGEN5016_S09	ZrGEN5016_S13					6		ARZXM5016N.MTN
	ZrGEN5018_S09	ZrGEN5018_S13					8		ARZXM5018N.MTN
	ZrGEN5025_S09	ZrGEN5025_S13				1.5	4.5		ARZXM5025N.MTN
	ZrGEN5026_S09	ZrGEN5026_S13					6		ARZXM5026N.MTN
	ZrGEN5028_S09	ZrGEN5028_S13					8		ARZXM5028N.MTN
	ZrGEN5035_S09	ZrGEN5035_S13				3.0	4.5		ARZXM5035N.MTN
	ZrGEN5036_S09	ZrGEN5036_S13					6		ARZXM5036N.MTN
	ZrGEN5038_S09	ZrGEN5038_S13					8		ARZXM5038N.MTN
	ZrGEN5045_S09	ZrGEN5045_S13				4.0	4.5		ARZXM5045N.MTN
	ZrGEN5046_S09	ZrGEN5046_S13					6		ARZXM5046N.MTN
	ZrGEN5048_S09	ZrGEN5048_S13					8		ARZXM5048N.MTN
	ZrGEN5515_S09	ZrGEN5515_S13	Core 4.8	Ø5.5	0.6	4.5	Hex	ARZXM5515.MTN	
	ZrGEN5516_S09	ZrGEN5516_S13				6		ARZXM5516.MTN	
	ZrGEN5518_S09	ZrGEN5518_S13				8		ARZXM5518.MTN	
	ZrGEN5525_S09	ZrGEN5525_S13			1.5	4.5		ARZXM5525.MTN	
	ZrGEN5526_S09	ZrGEN5526_S13				6		ARZXM5526.MTN	
	ZrGEN5528_S09	ZrGEN5528_S13				8		ARZXM5528.MTN	
	ZrGEN5535_S09	ZrGEN5535_S13			3.0	4.5		ARZXM5535.MTN	
	ZrGEN5536_S09	ZrGEN5536_S13				6		ARZXM5536.MTN	
	ZrGEN5538_S09	ZrGEN5538_S13				8		ARZXM5538.MTN	
	ZrGEN5545_S09	ZrGEN5545_S13			4.0	4.5		ARZXM5545.MTN	
	ZrGEN5546_S09	ZrGEN5546_S13				6		ARZXM5546.MTN	
	ZrGEN5548_S09	ZrGEN5548_S13				8		ARZXM5548.MTN	
	ZrGEN5515_S09	ZrGEN5515_S13		Core 4.8	Ø5.5	0.6	4.5	Non-hex	ARZXM5515N.MTN
	ZrGEN5516_S09	ZrGEN5516_S13					6		ARZXM5516N.MTN
	ZrGEN5518_S09	ZrGEN5518_S13					8		ARZXM5518N.MTN
	ZrGEN5525_S09	ZrGEN5525_S13				1.5	4.5		ARZXM5525N.MTN
	ZrGEN5526_S09	ZrGEN5526_S13					6		ARZXM5526N.MTN
	ZrGEN5528_S09	ZrGEN5528_S13					8		ARZXM5528N.MTN
	ZrGEN5535_S09	ZrGEN5535_S13				3.0	4.5		ARZXM5535N.MTN
	ZrGEN5536_S09	ZrGEN5536_S13					6		ARZXM5536N.MTN
	ZrGEN5538_S09	ZrGEN5538_S13					8		ARZXM5538N.MTN
	ZrGEN5545_S09	ZrGEN5545_S13				4.0	4.5		ARZXM5545N.MTN
	ZrGEN5546_S09	ZrGEN5546_S13					6		ARZXM5546N.MTN
	ZrGEN5548_S09	ZrGEN5548_S13					8		ARZXM5548N.MTN

AnyRidge\_Extra\_  
ZrGEN\_  
Custom abutment





**ZrGEN**

Ti-Base



Library info			Product info						
ZrGEN Option In 3shape library			Core (Ø)	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C	
System	Kit for 9mm	Kit for 13mm							
AnyRidge_Extra_ZrGEN_Custom abutment	ZrGEN5515_S09	ZrGEN5515_S13	Core 4.0	Ø5.0	0.6	4.5	Hex	ARZXL5515.MTN	
	ZrGEN5516_S09	ZrGEN5516_S13				6		ARZXL5516.MTN	
	ZrGEN5518_S09	ZrGEN5518_S13				8		ARZXL5518.MTN	
	ZrGEN5525_S09	ZrGEN5525_S13			1.5	4.5		ARZXL5525.MTN	
	ZrGEN5526_S09	ZrGEN5526_S13				6		ARZXL5526.MTN	
	ZrGEN5528_S09	ZrGEN5528_S13				8		ARZXL5528.MTN	
	ZrGEN5535_S09	ZrGEN5535_S13			3.0	4.5		ARZXL5535.MTN	
	ZrGEN5536_S09	ZrGEN5536_S13				6		ARZXL5536.MTN	
	ZrGEN5538_S09	ZrGEN5538_S13				8		ARZXL5538.MTN	
	ZrGEN5545_S09	ZrGEN5545_S13			4.0	4.5		ARZXL5545.MTN	
	ZrGEN5546_S09	ZrGEN5546_S13				6		ARZXL5546.MTN	
	ZrGEN5548_S09	ZrGEN5548_S13				8		ARZXL5548.MTN	
	ZrGEN5515_S09	ZrGEN5515_S13		Ø5.0	Core 4.0	0.6	4.5	Non-hex	ARZXL5515N.MTN
	ZrGEN5516_S09	ZrGEN5516_S13					6		ARZXL5516N.MTN
	ZrGEN5518_S09	ZrGEN5518_S13					8		ARZXL5518N.MTN
	ZrGEN5525_S09	ZrGEN5525_S13				1.5	4.5		ARZXL5525N.MTN
	ZrGEN5526_S09	ZrGEN5526_S13					6		ARZXL5526N.MTN
	ZrGEN5528_S09	ZrGEN5528_S13					8		ARZXL5528N.MTN
	ZrGEN5535_S09	ZrGEN5535_S13				3.0	4.5		ARZXL5535N.MTN
	ZrGEN5536_S09	ZrGEN5536_S13					6		ARZXL5536N.MTN
	ZrGEN5538_S09	ZrGEN5538_S13					8		ARZXL5538N.MTN
	ZrGEN5545_S09	ZrGEN5545_S13				4.0	4.5		ARZXL5545N.MTN
	ZrGEN5546_S09	ZrGEN5546_S13					6		ARZXL5546N.MTN
	ZrGEN5548_S09	ZrGEN5548_S13					8		ARZXL5548N.MTN
ZrGEN6015_S09	ZrGEN6015_S13	Ø5.5	Core 4.0	0.6	4.5	Hex	ARZXL6015.MTN		
ZrGEN6016_S09	ZrGEN6016_S13				6		ARZXL6016.MTN		
ZrGEN6018_S09	ZrGEN6018_S13				8		ARZXL6018.MTN		
ZrGEN6025_S09	ZrGEN6025_S13			1.5	4.5		ARZXL6025.MTN		
ZrGEN6026_S09	ZrGEN6026_S13				6		ARZXL6026.MTN		
ZrGEN6028_S09	ZrGEN6028_S13				8		ARZXL6028.MTN		
ZrGEN6035_S09	ZrGEN6035_S13			3.0	4.5		ARZXL6035.MTN		
ZrGEN6036_S09	ZrGEN6036_S13				6		ARZXL6036.MTN		
ZrGEN6038_S09	ZrGEN6038_S13				8		ARZXL6038.MTN		
ZrGEN6045_S09	ZrGEN6045_S13			4.0	4.5		ARZXL6045.MTN		
ZrGEN6046_S09	ZrGEN6046_S13				6		ARZXL6046.MTN		
ZrGEN6048_S09	ZrGEN6048_S13				8		ARZXL6048.MTN		
ZrGEN6015_S09	ZrGEN6015_S13	Ø5.5	Core 4.0	0.6	4.5	Non-hex	ARZXL6015N.MTN		
ZrGEN6016_S09	ZrGEN6016_S13				6		ARZXL6016N.MTN		
ZrGEN6018_S09	ZrGEN6018_S13				8		ARZXL6018N.MTN		
ZrGEN6025_S09	ZrGEN6025_S13			1.5	4.5		ARZXL6025N.MTN		
ZrGEN6026_S09	ZrGEN6026_S13				6		ARZXL6026N.MTN		
ZrGEN6028_S09	ZrGEN6028_S13				8		ARZXL6028N.MTN		
ZrGEN6035_S09	ZrGEN6035_S13			3.0	4.5		ARZXL6035N.MTN		
ZrGEN6036_S09	ZrGEN6036_S13				6		ARZXL6036N.MTN		
ZrGEN6038_S09	ZrGEN6038_S13				8		ARZXL6038N.MTN		
ZrGEN6045_S09	ZrGEN6045_S13			4.0	4.5		ARZXL6045N.MTN		
ZrGEN6046_S09	ZrGEN6046_S13				6		ARZXL6046N.MTN		
ZrGEN6048_S09	ZrGEN6048_S13				8		ARZXL6048N.MTN		





Pre-milled



Bevel Platform




AnyRidge\_Extra\_TiGEN  
\_Custom abutment

Library info			Product info				
Pre-milled Option In 3shape library			Core (∅)	Diameter (∅)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
AnyRidge_Extra_TiGEN _Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	3.3 Bevel	10	20	Hex	ARTXN1020.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		ARTXN1020N.MTN
	TiGEN 1020_S09	TiGEN 1020_S13		10	20	Non-Hex	ARTXN1220.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		ARTXN1220N.MTN
	TiGEN 1020_S09	TiGEN 1020_S13	4.0 Platform	10	20	Hex	ARTXM1020.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		ARTXM1020N.MTN
	TiGEN 1020_S09	TiGEN 1020_S13		10	20	Non-Hex	ARTXM1220.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		ARTXM1220N.MTN
	TiGEN 1020_S09	TiGEN 1020_S13	4.8 Platform	10	20	Hex	ARTXL1020.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		ARTXL1020N.MTN
	TiGEN 1020_S09	TiGEN 1020_S13		10	20	Non-Hex	ARTXL1220.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		ARTXL1220N.MTN




# AnyOne



\* EZ Post Prosthesis by utilizing the stock abutment like a scan abutment.  
\*\* C-type For Sirona CEREC System


# Scan abutment

## Scan Abutment [AO]



Product info			Library info (Choose below CAD library option)		
Diameter	Lenght	Ref.C	exocad library name	3shape library name	Dental wings library name
Ø4.0	9	AAOISR4009T	AnyOne Custom Abutment (9mm)	AnyOne_ZrGEN_Custom abutment	AnyOne Internal
	13	AAOISR4013T	AnyOne Custom Abutment (13mm)	AnyOne_TiGEN_Custom abutment	

## Scan Healing Abutment [AO]



Product info				Library info (Choose below CAD library option)		
Diameter (D)	Height	Abutment Ref.C	Scan Post Ref.C	exocad library name	3shape library name	Dental wings library name
Ø4.0	4	AOISH4004T	SP4007.MTN	AnyOne SHA Custom Abutment (Healing Height 4mm)	For ZrGEN AnyOne_ZrGEN_SHA_Custom abutment	
	5	AOISH4005T		AnyOne SHA Custom Abutment (Healing Height 5mm)		
	7	AOISH4007T		AnyOne SHA Custom Abutment (Healing Height 7mm)		
Ø4.5	4	AOISH4504T	SP5007.MTN	AnyOne SHA Custom Abutment (Healing Height 4mm)	ZrGEN4015_Healing Height 4mm ZrGEN4015_Healing Height 5mm ZrGEN4015_Healing Height 7mm	
	5	AOISH4505T		AnyOne SHA Custom Abutment (Healing Height 5mm)		
	7	AOISH4507T		AnyOne SHA Custom Abutment (Healing Height 7mm)		
Ø5.5	4	AOISH5504T	SP6007.MTN	AnyOne SHA Custom Abutment (Healing Height 4mm)	For TiGEN AnyOne_TiGEN_SHA_Custom abutment	
	5	AOISH5505T		AnyOne SHA Custom Abutment (Healing Height 5mm)		
	7	AOISH5507T		AnyOne SHA Custom Abutment (Healing Height 7mm)		
Ø6.5	4	AOISH6504T	SP7007.MTN	AnyOne SHA Custom Abutment (Healing Height 4mm)	TiGEN1020_Healing Height 4mm TiGEN1020_Healing Height 5mm TiGEN1020_Healing Height 7mm	
	5	AOISH6505T		AnyOne SHA Custom Abutment (Healing Height 5mm)		
	7	AOISH6507T		AnyOne SHA Custom Abutment (Healing Height 7mm)		

## Scan Abutment C-type [AO] For use CEREC



Product info				
Diameter	Cuff Height	Post height	Post size	Ref. C
Ø3.9	0.5	4.7	Small	AOICSS3405T
	1.0			AOICSS3410T
	2.0			AOICSS3420T
Ø4.3	0.5			AOICSS3805T
	1.0			AOICSS3810T
	2.0			AOICSS3820T
Ø5.5	0.5		Large	AOICSL4505T
	1.0			AOICSL4510T
	2.0			AOICSL4520T

\* Compatible with [Xive Library](#) in **CEREC Or in-Lab** CAD Software.

**ZrGEN**  
Ti-Base



Library info		Product info							
ZrGEN Option In exocad library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C			
System	Kit								
ZrGEN	ZrGEN 4015	Ø4.0	0.6	4.5	Hex	AAOIPR4015.MTN			
	ZrGEN 4016			6.0		AAOIPR4016.MTN			
	ZrGEN 4018			8.0		AAOIPR4018.MTN			
	ZrGEN 4025		1.5	4.5		AAOIPR4025.MTN			
	ZrGEN 4026			6.0		AAOIPR4026.MTN			
	ZrGEN 4028			8.0		AAOIPR4028.MTN			
	ZrGEN 4035		3.0	4.5		AAOIPR4035.MTN			
	ZrGEN 4036			6.0		AAOIPR4036.MTN			
	ZrGEN 4038			8.0		AAOIPR4038.MTN			
	ZrGEN 4045		4.0	4.5		AAOIPR4045.MTN			
	ZrGEN 4046			6.0		AAOIPR4046.MTN			
	ZrGEN 4048			8.0		AAOIPR4048.MTN			
	AnyOne Custom Abutment (9, 13mm)		ZrGEN 4015	Ø4.0		0.6	4.5	Non-Hex	AAOIPR4015N.MTN
			ZrGEN 4016				6.0		AAOIPR4016N.MTN
ZrGEN 4018		8.0	AAOIPR4018N.MTN						
ZrGEN 4025		1.5	4.5		AAOIPR4025N.MTN				
ZrGEN 4026			6.0		AAOIPR4026N.MTN				
ZrGEN 4028			8.0		AAOIPR4028N.MTN				
ZrGEN 4035		3.0	4.5		AAOIPR4035N.MTN				
ZrGEN 4036			6.0		AAOIPR4036N.MTN				
ZrGEN 4038			8.0		AAOIPR4038N.MTN				
ZrGEN 4045		4.0	4.5		AAOIPR4045N.MTN				
ZrGEN 4046			6.0		AAOIPR4046N.MTN				
ZrGEN 4048			8.0		AAOIPR4048N.MTN				
AnyOne SHA Custom Abutment (Healing Height 4, 5, 7mm)		ZrGEN 4515	Ø4.5		0.6	4.5	Hex		AAOIPR4515.MTN
		ZrGEN 4516				6.0			AAOIPR4516.MTN
	ZrGEN 4518	8.0		AAOIPR4518.MTN					
	ZrGEN 4525	1.5		4.5	AAOIPR4525.MTN				
	ZrGEN 4526			6.0	AAOIPR4526.MTN				
	ZrGEN 4528			8.0	AAOIPR4528.MTN				
	ZrGEN 4535	3.0		4.5	AAOIPR4535.MTN				
	ZrGEN 4536			6.0	AAOIPR4536.MTN				
	ZrGEN 4538			8.0	AAOIPR4538.MTN				
	ZrGEN 4545	4.0		4.5	AAOIPR4545.MTN				
	ZrGEN 4546			6.0	AAOIPR4546.MTN				
	ZrGEN 4548			8.0	AAOIPR4548.MTN				
	AnyOne SHA Custom Abutment (Healing Height 4, 5, 7mm)	ZrGEN 4515		Ø4.5	0.6	4.5		Non-Hex	AAOIPR4515N.MTN
		ZrGEN 4516				6.0			AAOIPR4516N.MTN
ZrGEN 4518		8.0	AAOIPR4518N.MTN						
ZrGEN 4525		1.5	4.5		AAOIPR4525N.MTN				
ZrGEN 4526			6.0		AAOIPR4526N.MTN				
ZrGEN 4528			8.0		AAOIPR4528N.MTN				
ZrGEN 4535		3.0	4.5		AAOIPR4535N.MTN				
ZrGEN 4536			6.0		AAOIPR4536N.MTN				
ZrGEN 4538			8.0		AAOIPR4538N.MTN				
ZrGEN 4545		4.0	4.5		AAOIPR4545N.MTN				
ZrGEN 4546			6.0		AAOIPR4546N.MTN				
ZrGEN 4548			8.0		AAOIPR4548N.MTN				

\* ZrGEN Library provides Cement gap (ZrGEN and Abutment space) 0, 10, 20, 30, 40, 50 µ.



**TiGEN**  
Pre-milled



Library info		Product info			
Pre-milled Option In exocad library		Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit				
<p>AnyOne Custom Abutment (9, 13mm)</p> <p>AnyOne SHA Custom Abutment (Healing Height 4, 5, 7mm)</p>	TiGEN 1020	10	20	Hex	AOTR1020.MTN
	TiGEN 1220	12	20		AOTR1220.MTN
	TiGEN 1020	10	20	Non-Hex	AOTR1020N.MTN
	TiGEN 1220	12	20		AOTR1220N.MTN



# EZ Post

Stock abutment





Library info		Product info						
EZ Post Option In exocad library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C		
System	Kit							
	EZ Post 4511	Ø4.5	1.0	4.0	Hex	EP4511HT		
	EZ Post 4514		1.5			EP4514HT		
	EZ Post 4524		2.5			EP4524HT		
	EZ Post 4534		3.5	EP4534HT				
	EZ Post 4544		4.5	EP4544HT				
	EZ Post 4554		5.5	EP4554HT				
	EZ Post 4510		1.0	5.5		EP4510HT		
	EZ Post 4515		1.5			EP4515HT		
	EZ Post 4525		2.5			EP4525HT		
	EZ Post 4535		3.5	EP4535HT				
	EZ Post 4545		4.5	EP4545HT				
	EZ Post 4555		5.5	EP4555HT				
	EZ Post 4517		1.5	7.0		EP4517HT		
	EZ Post 4527		2.5			EP4527HT		
	EZ Post 4537		3.5			EP4537HT		
	EZ Post 4547		4.5	EP4547HT				
	EZ Post 4557		5.5	EP4557HT				
	EZ Post 4511	Ø4.5	1.0	4.0	Non-Hex	EP4511NT		
	EZ Post 4514		1.5			EP4514NT		
	EZ Post 4524		2.5			EP4524NT		
	EZ Post 4534		3.5	EP4534NT				
	EZ Post 4544		4.5	EP4544NT				
	EZ Post 4554		5.5	EP4554NT				
	EZ Post 4510		1.0	5.0		EP4510NT		
	EZ Post 4515		1.5			EP4515NT		
	EZ Post 4525		2.5			EP4525NT		
	EZ Post 4535		3.5	EP4535NT				
	EZ Post 4545		4.5	EP4545NT				
	EZ Post 4555		5.5	EP4555NT				
	EZ Post 4517		1.5	7.0		EP4517NT		
	EZ Post 4527		2.5			EP4527NT		
	EZ Post 4537		3.5			EP4537NT		
	EZ Post 4547		4.5	EP4547NT				
	EZ Post 4557		5.5	EP4557NT				
	EZ Post 5514	Ø5.5	1.5	4.0	Hex	EP5514HT		
	EZ Post 5524		2.5			EP5524HT		
	EZ Post 5534		3.5			EP5534HT		
	EZ Post 5544		4.5	EP5544HT				
	EZ Post 5554		5.5	EP5554HT				
	EZ Post 5515		1.5	5.5		EP5515HT		
	EZ Post 5525		2.5			EP5525HT		
	EZ Post 5535		3.5			EP5535HT		
	EZ Post 5545		4.5	EP5545HT				
	EZ Post 5555		5.5	EP5555HT				
	EZ Post 5517		1.5	7.0		EP5517HT		
	EZ Post 5527		2.5			EP5527HT		
	EZ Post 5537		3.5			EP5537HT		
	EZ Post 5547		4.5	EP5547HT				
	EZ Post 5557		5.5	EP5557HT				
	EZ Post 5514		Ø5.5	1.5		4.0	Non-Hex	EP5514NT
	EZ Post 5524			2.5				EP5524NT
	EZ Post 5534	3.5		EP5534NT				
	EZ Post 5544	4.5		EP5544NT				
	EZ Post 5554	5.5		EP5554NT				
	EZ Post 5515	1.5		5.5	EP5515NT			
	EZ Post 5525	2.5			EP5525NT			
	EZ Post 5535	3.5			EP5535NT			
	EZ Post 5545	4.5		EP5545NT				
	EZ Post 5555	5.5		EP5555NT				
	EZ Post 5517	1.5		7.0	EP5517NT			
	EZ Post 5527	2.5			EP5527NT			
	EZ Post 5537	3.5			EP5537NT			
	EZ Post 5547	4.5		EP5547NT				
	EZ Post 5557	5.5		EP5557NT				
	EZ Post 6514	Ø6.5		1.5	4.0	Hex		EP6514HT
	EZ Post 6524			2.5				EP6524HT
	EZ Post 6534		3.5	EP6534HT				
	EZ Post 6544		4.5	EP6544HT				
	EZ Post 6554		5.5	EP6554HT				
	EZ Post 6515		1.5	5.5	EP6515HT			
	EZ Post 6525		2.5		EP6525HT			
	EZ Post 6535		3.5		EP6535HT			
	EZ Post 6545		4.5	EP6545HT				
	EZ Post 6555		5.5	EP6555HT				
	EZ Post 6157		1.5	7.0	EP6157HT			
	EZ Post 6527		2.5		EP6527HT			
	EZ Post 6537		3.5		EP6537HT			
	EZ Post 6547		4.5	EP6547HT				
	EZ Post 6557		5.5	EP6557HT				
	EZ Post 6514		Ø6.5	1.5	4.0		Non-Hex	EP6514NT
	EZ Post 6524			2.5				EP6524NT
	EZ Post 6534	3.5		EP6534NT				
	EZ Post 6544	4.5		EP6544NT				
	EZ Post 6554	5.5		EP6554NT				
	EZ Post 6515	1.5		5.5	EP6515NT			
	EZ Post 6525	2.5			EP6525NT			
	EZ Post 6535	3.5			EP6535NT			
	EZ Post 6545	4.5		EP6545NT				
	EZ Post 6555	5.5		EP6555NT				
	EZ Post 6157	1.5		7.0	EP6157NT			
	EZ Post 6527	2.5			EP6527NT			
	EZ Post 6537	3.5			EP6537NT			
	EZ Post 6547	4.5		EP6547NT				
	EZ Post 6557	5.5		EP6557NT				



**ZrGEN**

Ti-Base



Library info			Product info							
ZrGEN Option In 3shape library			Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C			
System	Kit for 9mm	Kit for 13mm								
	ZrGEN 4015_S09	ZrGEN 4015_S13	Ø4.0	0.6	4.5	Hex	AAOIPR4015.MTN			
	ZrGEN 4016_S09	ZrGEN 4016_S13			6.0		AAOIPR4016.MTN			
	ZrGEN 4018_S09	ZrGEN 4018_S13			8.0		AAOIPR4018.MTN			
	ZrGEN 4025_S09	ZrGEN 4025_S13			1.5		4.5	AAOIPR4025.MTN		
	ZrGEN 4026_S09	ZrGEN 4026_S13					6.0	AAOIPR4026.MTN		
	ZrGEN 4028_S09	ZrGEN 4028_S13					8.0	AAOIPR4028.MTN		
	ZrGEN 4035_S09	ZrGEN 4035_S13			3.0		4.5	AAOIPR4035.MTN		
	ZrGEN 4036_S09	ZrGEN 4036_S13					6.0	AAOIPR4036.MTN		
	ZrGEN 4038_S09	ZrGEN 4038_S13					8.0	AAOIPR4038.MTN		
	ZrGEN 4045_S09	ZrGEN 4045_S13		4.0	4.5		AAOIPR4045.MTN			
	ZrGEN 4046_S09	ZrGEN 4046_S13			6.0		AAOIPR4046.MTN			
	ZrGEN 4048_S09	ZrGEN 4048_S13			8.0		AAOIPR4048.MTN			
	ZrGEN 4015_S09	ZrGEN 4015_S13		Ø4.0	0.6		4.5	Non-Hex	AAOIPR4015N.MTN	
	ZrGEN 4016_S09	ZrGEN 4016_S13					6.0		AAOIPR4016N.MTN	
	ZrGEN 4018_S09	ZrGEN 4018_S13					8.0		AAOIPR4018N.MTN	
	ZrGEN 4025_S09	ZrGEN 4025_S13					1.5		4.5	AAOIPR4025N.MTN
	ZrGEN 4026_S09	ZrGEN 4026_S13							6.0	AAOIPR4026N.MTN
	ZrGEN 4028_S09	ZrGEN 4028_S13							8.0	AAOIPR4028N.MTN
	ZrGEN 4035_S09	ZrGEN 4035_S13					3.0		4.5	AAOIPR4035N.MTN
	ZrGEN 4036_S09	ZrGEN 4036_S13							6.0	AAOIPR4036N.MTN
	ZrGEN 4038_S09	ZrGEN 4038_S13							8.0	AAOIPR4038N.MTN
	ZrGEN 4045_S09	ZrGEN 4045_S13			4.0		4.5		AAOIPR4045N.MTN	
	ZrGEN 4046_S09	ZrGEN 4046_S13					6.0		AAOIPR4046N.MTN	
	ZrGEN 4048_S09	ZrGEN 4048_S13					8.0		AAOIPR4048N.MTN	
	ZrGEN 4515_S09	ZrGEN 4515_S13	Ø4.5		0.6	4.5	Hex		AAOIPR4515.MTN	
	ZrGEN 4516_S09	ZrGEN 4516_S13				6.0			AAOIPR4516.MTN	
	ZrGEN 4518_S09	ZrGEN 4518_S13				8.0			AAOIPR4518.MTN	
	ZrGEN 4525_S09	ZrGEN 4525_S13				1.5			4.5	AAOIPR4525.MTN
	ZrGEN 4526_S09	ZrGEN 4526_S13							6.0	AAOIPR4526.MTN
	ZrGEN 4528_S09	ZrGEN 4528_S13							8.0	AAOIPR4528.MTN
	ZrGEN 4535_S09	ZrGEN 4535_S13				3.0			4.5	AAOIPR4535.MTN
	ZrGEN 4536_S09	ZrGEN 4536_S13							6.0	AAOIPR4536.MTN
	ZrGEN 4538_S09	ZrGEN 4538_S13							8.0	AAOIPR4538.MTN
	ZrGEN 4545_S09	ZrGEN 4545_S13			4.0	4.5			AAOIPR4545.MTN	
	ZrGEN 4546_S09	ZrGEN 4546_S13				6.0			AAOIPR4546.MTN	
	ZrGEN 4548_S09	ZrGEN 4548_S13				8.0			AAOIPR4548.MTN	
	ZrGEN 4515_S09	ZrGEN 4515_S13		Ø4.5	0.6	4.5		Non-Hex	AAOIPR4515N.MTN	
	ZrGEN 4516_S09	ZrGEN 4516_S13				6.0			AAOIPR4516N.MTN	
	ZrGEN 4518_S09	ZrGEN 4518_S13				8.0			AAOIPR4518N.MTN	
	ZrGEN 4525_S09	ZrGEN 4525_S13				1.5			4.5	AAOIPR4525N.MTN
	ZrGEN 4526_S09	ZrGEN 4526_S13							6.0	AAOIPR4526N.MTN
	ZrGEN 4528_S09	ZrGEN 4528_S13							8.0	AAOIPR4528N.MTN
	ZrGEN 4535_S09	ZrGEN 4535_S13				3.0			4.5	AAOIPR4535N.MTN
	ZrGEN 4536_S09	ZrGEN 4536_S13							6.0	AAOIPR4536N.MTN
	ZrGEN 4538_S09	ZrGEN 4538_S13							8.0	AAOIPR4538N.MTN
	ZrGEN 4545_S09	ZrGEN 4545_S13			4.0	4.5			AAOIPR4545N.MTN	
	ZrGEN 4546_S09	ZrGEN 4546_S13				6.0			AAOIPR4546N.MTN	
	ZrGEN 4548_S09	ZrGEN 4548_S13				8.0			AAOIPR4548N.MTN	

\* ZrGEN Library provides Cement gap (ZrGEN with Zirconia abutment space) 0, 10, 20, 30, 40, 50 µ.





**TiGEN**  
Pre-milled



Library info			Product info			
Pre-milled Option In 3shape library			Diameter (∅)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm				
 AnyOne_TiGEN_ Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	10	20	Hex	AOTR1020.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	12	20		AOTR1220.MTN
 AnyOne_TiGEN_SHA Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	10	20	Non-Hex	AOTR1020N.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	12	20		AOTR1220N.MTN



# EZ Post


Stock abutment




Library info		Product info						
EZ Post Option In 3shape library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C		
System	Kit							
	EZ Post 4511	Ø4.5	1.0	4.0	Hex	EP4511HT		
	EZ Post 4514		1.5			EP4514HT		
	EZ Post 4524		2.5			EP4524HT		
	EZ Post 4534		3.5	EP4534HT				
	EZ Post 4544		4.5	EP4544HT				
	EZ Post 4554		5.5	EP4554HT				
	EZ Post 4510		1.0	5.5		EP4510HT		
	EZ Post 4515		1.5			EP4515HT		
	EZ Post 4525		2.5			EP4525HT		
	EZ Post 4535		3.5	EP4535HT				
	EZ Post 4545		4.5	EP4545HT				
	EZ Post 4555		5.5	EP4555HT				
	EZ Post 4517		1.5	7.0		EP4517HT		
	EZ Post 4527		2.5			EP4527HT		
	EZ Post 4537		3.5			EP4537HT		
	EZ Post 4547		4.5	EP4547HT				
	EZ Post 4557		5.5	EP4557HT				
	EZ Post 4511	Ø4.5	1.0	4.0	Non-Hex	EP4511NT		
	EZ Post 4514		1.5			EP4514NT		
	EZ Post 4524		2.5			EP4524NT		
	EZ Post 4534		3.5	EP4534NT				
	EZ Post 4544		4.5	EP4544NT				
	EZ Post 4554		5.5	EP4554NT				
	EZ Post 4510		1.0	5.0		EP4510NT		
	EZ Post 4515		1.5			EP4515NT		
	EZ Post 4525		2.5			EP4525NT		
	EZ Post 4535		3.5	EP4535NT				
	EZ Post 4545		4.5	EP4545NT				
	EZ Post 4555		5.5	EP4555NT				
	EZ Post 4517		1.5	7.0		EP4517NT		
	EZ Post 4527		2.5			EP4527NT		
	EZ Post 4537		3.5			EP4537NT		
	EZ Post 4547		4.5	EP4547NT				
	EZ Post 4557		5.5	EP4557NT				
	EZ Post 5514	Ø5.5	1.5	4.0	Hex	EP5514HT		
	EZ Post 5524		2.5			EP5524HT		
	EZ Post 5534		3.5			EP5534HT		
	EZ Post 5544		4.5	EP5544HT				
	EZ Post 5554		5.5	EP5554HT				
	EZ Post 5515		1.5	5.5		EP5515HT		
	EZ Post 5525		2.5			EP5525HT		
	EZ Post 5535		3.5			EP5535HT		
	EZ Post 5545		4.5	EP5545HT				
	EZ Post 5555		5.5	EP5555HT				
	EZ Post 5517		1.5	7.0		EP5517HT		
	EZ Post 5527		2.5			EP5527HT		
	EZ Post 5537		3.5			EP5537HT		
	EZ Post 5547		4.5	EP5547HT				
	EZ Post 5557		5.5	EP5557HT				
	EZ Post 5514		Ø5.5	1.5		4.0	Non-Hex	EP5514NT
	EZ Post 5524			2.5				EP5524NT
	EZ Post 5534	3.5		EP5534NT				
	EZ Post 5544	4.5		EP5544NT				
	EZ Post 5554	5.5		EP5554NT				
	EZ Post 5515	1.5		5.5	EP5515NT			
	EZ Post 5525	2.5			EP5525NT			
	EZ Post 5535	3.5			EP5535NT			
	EZ Post 5545	4.5		EP5545NT				
	EZ Post 5555	5.5		EP5555NT				
	EZ Post 5517	1.5		7.0	EP5517NT			
	EZ Post 5527	2.5			EP5527NT			
	EZ Post 5537	3.5			EP5537NT			
	EZ Post 5547	4.5		EP5547NT				
	EZ Post 5557	5.5		EP5557NT				
	EZ Post 6514	Ø6.5		1.5	4.0	Hex		EP6514HT
	EZ Post 6524			2.5				EP6524HT
	EZ Post 6534		3.5	EP6534HT				
	EZ Post 6544		4.5	EP6544HT				
	EZ Post 6554		5.5	EP6554HT				
	EZ Post 6515		1.5	5.5	EP6515HT			
	EZ Post 6525		2.5		EP6525HT			
	EZ Post 6535		3.5		EP6535HT			
	EZ Post 6545		4.5	EP6545HT				
	EZ Post 6555		5.5	EP6555HT				
	EZ Post 6157		1.5	7.0	EP6157HT			
	EZ Post 6527		2.5		EP6527HT			
	EZ Post 6537		3.5		EP6537HT			
	EZ Post 6547		4.5	EP6547HT				
	EZ Post 6557		5.5	EP6557HT				
	EZ Post 6514		Ø6.5	1.5	4.0		Non-Hex	EP6514NT
	EZ Post 6524			2.5				EP6524NT
	EZ Post 6534	3.5		EP6534NT				
	EZ Post 6544	4.5		EP6544NT				
	EZ Post 6554	5.5		EP6554NT				
	EZ Post 6515	1.5		5.5	EP6515NT			
	EZ Post 6525	2.5			EP6525NT			
	EZ Post 6535	3.5			EP6535NT			
	EZ Post 6545	4.5		EP6545NT				
	EZ Post 6555	5.5		EP6555NT				
	EZ Post 6157	1.5		7.0	EP6157NT			
	EZ Post 6527	2.5			EP6527NT			
	EZ Post 6537	3.5			EP6537NT			
	EZ Post 6547	4.5		EP6547NT				
	EZ Post 6557	5.5		EP6557NT				



AnyOne\_EZ Post\_Scan abutment



AnyOne\_EZ Post\_SHA




AnOne\_EZ Post\_EZ Post



**ZrGEN**

Ti-Base




Library info			Product info				
ZrGEN Option In Dental wings library			Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
 AnyOne Internal ZrGEN (A_07A6)	ZrGEN 4015_S09	ZrGEN 4015_S13	Ø4.0	0.6	4.5	Hex	AAOIPR4015.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13			6		AAOIPR4016.MTN
	ZrGEN 4015_S09	ZrGEN 4015_S13	Ø4.0	0.6	4.5	Non-Hex	AAOIPR4015N.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13			6		AAOIPR4016N.MTN
	ZrGEN 4525_S09	ZrGEN 4525_S13	Ø4.5	0.6	4.5	Hex	AAOIPR4515.MTN
	ZrGEN 4526_S09	ZrGEN 4526_S13			6		AAOIPR4516.MTN
	ZrGEN 4525_S09	ZrGEN 4525_S13	Ø4.5	0.6	4.5	Non-Hex	AAOIPR4515N.MTN
	ZrGEN 4526_S09	ZrGEN 4526_S13			6		AAOIPR4516N.MTN



**TiGEN**  
Pre-milled



Library info			Product info			
Pre-milled Option In Dental wings library			Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm				
 AnyOne Internal ↓ TiGEN (A_07A6)	TiGEN 1220_S09	TiGEN 1220_S13	12	20	Hex	AOTR1220.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	12	20	Non-Hex	AOTR1220N.MTN




**ZrGEN**  
For CEREC



Product info				
Diameter (Ø)	Cuff height (mm)	Post height (mm)	Post size (mm)	Ref. C
Ø3.9	0.5	4.7	Small	AOCS3405.MTN
	1			AOCS3410.MTN
	2			AOCS3420.MTN
Ø4.3	0.5			AOCS3805.MTN
	1			AOCS3810.MTN
	2			AOCS3820.MTN
Ø5.5	0.5		Large	AOCL4505.MTN
	1			AOCL4510.MTN
	2			AOCL4520.MTN


\* Compatible with [Xive Library](#) in **CEREC Or in-Lab** CAD Software.

# AnyOne OneStage



\* **EZ Post** We can use the EZ Post as a Scan abutment. For user convenience.

# Scan abutment



**Scan Abutment [AOO]** \* Compatible with [EXi\(Exfeel internal\)](#) Scan abutment

Product info			Library info (Choose below CAD library option)		
Diameter	Lenght	Ref.C	exocad library name	3shape library name	Dental wings library name
Ø4.0	10	AEXISR4010T	AnyOne Onestage Custom Abutment (10mm)	AnyOne Onestage_ZrGEN_Custom abutment AnyOne Onestage_TIGEN_Custom abutment	-



**ZrGEN**

Ti-Base



AnyOne OneStage  
Custom Abutment (10mm)

Library info		Product info				
ZrGEN Option In exocad library		Diameter (∅)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
AnyOne OneStage Custom Abutment (10mm)	ZrGEN 5015	∅4.8	0.6	4.5	Octa	AEXIPR5015.MTN
	ZrGEN 5016			6.0		AEXIPR5016.MTN
	ZrGEN 5018			8.0		AEXIPR5018.MTN
	ZrGEN 5025		1.5	4.5		AEXIPR5025.MTN
	ZrGEN 5026			6.0		AEXIPR5026.MTN
	ZrGEN 5028			8.0		AEXIPR5028.MTN
	ZrGEN 5035		3.0	4.5		AEXIPR5035.MTN
	ZrGEN 5036			6.0		AEXIPR5036.MTN
	ZrGEN 5038			8.0		AEXIPR5038.MTN
	ZrGEN 5045		4.0	4.5		AEXIPR5045.MTN
	ZrGEN 5046			6.0		AEXIPR5046.MTN
	ZrGEN 5048			8.0		AEXIPR5048.MTN



\* ZrGEN Library provides Cement gap (ZrGEN with Zirconia Abutment space) 0, 10, 20, 30, 40, 50 μ.





**EZ Post**  
Stock abutment




Library info		Product info				
EZ Post Option In exocad library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
 <p>AnyOne Onestage Scanbody - EZpost (10mm)</p>	EZ Post 5004	Ø4.0	0.0	4.0	Hex	AOOEPO5004T
	EZ Post 5014		1.0			AOOEPO5014T
	EZ Post 5024		2.0			AOOEPO5024T
	EZ Post 5034		3.0			AOOEPO5034T
	EZ Post 5005	Ø5.5	0.0	5.5		AOOEPO5005T
	EZ Post 5015		1.0			AOOEPO5015T
	EZ Post 5025		2.0			AOOEPO5025T
	EZ Post 5035		3.0			AOOEPO5035T
	EZ Post 5007	Ø7.0	0.0	7.0		AOOEPO5007T
	EZ Post 5017		1.0			AOOEPO5017T
	EZ Post 5027		2.0			AOOEPO5027T
	EZ Post 5037		3.0			AOOEPO5037T
 <p>AnyOne Onestage EZpost - EZpost</p>	EZ Post 5004	Ø4.0	0.0	4.0	Non-Hex	AOOEPO5004T
	EZ Post 5014		1.0			AOOEPO5014T
	EZ Post 5024		2.0			AOOEPO5024T
	EZ Post 5034		3.0			AOOEPO5034T
	EZ Post 5005	Ø5.5	0.0	5.5		AOOEPO5005T
	EZ Post 5015		1.0			AOOEPO5015T
	EZ Post 5025		2.0			AOOEPO5025T
	EZ Post 5035		3.0			AOOEPO5035T
	EZ Post 5007	Ø7.0	0.0	7.0		AOOEPO5007T
	EZ Post 5017		1.0			AOOEPO5017T
	EZ Post 5027		2.0			AOOEPO5027T
	EZ Post 5037		3.0			AOOEPO5037T



**ZrGEN**  
Ti-Base





Library info		Product info				
ZrGEN Option In 3shape library		Diameter (∅)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 10mm					
 <p>AnyOne_OneStage_ZrGEN_ Custom abutment</p>	ZrGEN 5015_S10	∅4.8	0.6	4.5	Octa	AEXIPR5015.MTN
	ZrGEN 5016_S10			6.0		AEXIPR5016.MTN
	ZrGEN 5018_S10			8.0		AEXIPR5018.MTN
	ZrGEN 5025_S10		1.5	4.5		AEXIPR5025.MTN
	ZrGEN 5026_S10			6.0		AEXIPR5026.MTN
	ZrGEN 5028_S10			8.0		AEXIPR5028.MTN
	ZrGEN 5035_S10		3.0	4.5		AEXIPR5035.MTN
	ZrGEN 5036_S10			6.0		AEXIPR5036.MTN
	ZrGEN 5038_S10			8.0		AEXIPR5038.MTN
	ZrGEN 5045_S10		4.0	4.5		AEXIPR5045.MTN
	ZrGEN 5046_S10			6.0		AEXIPR5046.MTN
	ZrGEN 5048_S10			8.0		AEXIPR5048.MTN




**EZ Post**  
Stock abutment



Library info		Product info				
EZ Post Option In 3shape library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
 <p>AnyOne Onestage Scanbody - EZpost (10mm)</p>	EZ Post 5004	Ø4.0	0.0	4.0	Hex	AOOEPO5004T
	EZ Post 5014		1.0			AOOEPO5014T
	EZ Post 5024		2.0			AOOEPO5024T
	EZ Post 5034		3.0			AOOEPO5034T
	EZ Post 5005	Ø5.5	0.0	5.5		AOOEPO5005T
	EZ Post 5015		1.0			AOOEPO5015T
	EZ Post 5025		2.0			AOOEPO5025T
	EZ Post 5035		3.0			AOOEPO5035T
	EZ Post 5007	Ø7.0	0.0	7.0		AOOEPO5007T
	EZ Post 5017		1.0			AOOEPO5017T
	EZ Post 5027		2.0			AOOEPO5027T
	EZ Post 5037		3.0			AOOEPO5037T
 <p>AnyOne Onestage EZpost - EZpost</p>	EZ Post 5004	Ø4.0	0.0	4.0	Non-Hex	AOOEPO5004T
	EZ Post 5014		1.0			AOOEPO5014T
	EZ Post 5024		2.0			AOOEPO5024T
	EZ Post 5034		3.0			AOOEPO5034T
	EZ Post 5005	Ø5.5	0.0	5.5		AOOEPO5005T
	EZ Post 5015		1.0			AOOEPO5015T
	EZ Post 5025		2.0			AOOEPO5025T
	EZ Post 5035		3.0			AOOEPO5035T
	EZ Post 5007	Ø7.0	0.0	7.0		AOOEPO5007T
	EZ Post 5017		1.0			AOOEPO5017T
	EZ Post 5027		2.0			AOOEPO5027T
	EZ Post 5037		3.0			AOOEPO5037T




# AnyOne External



# Scan abutment

Scan Abutment [AOE] \* Compatible **Regular** with **Wide** Scan abutment



Product info				Library info (Choose below CAD library option)	
Diameter	Lenght	Type	Ref.C	exocad library name	3shape library name
Ø4.0	9	Small	AEXESS4009T	AnyOne_External_Small_Custom Abutment (9mm)	AnyOne_Extennaal_Small ZrGEN_Custom abutment
	13		AEXESS4013T	AnyOne_External_Small_Custom Abutment (13mm)	AnyOne_Extennaal_Small TiGEN_Custom abutment
Ø4.0	9	Regular	AEXESR4009T	AnyOne_External_Regular_Custom Abutment (9mm)	AnyOne_Extennaal_Regular ZrGEN_Custom abutment
	13		AEXESR4013T	AnyOne_External_Regular_Custom Abutment (13mm)	AnyOne_Extennaal_Regular TiGEN_Custom abutment



**ZrGEN**  
Ti-Base



AnyOne External  
Small Custom Abutment  
(9, 13mm)



AnyOne External  
Regular Custom Abutment  
(9, 13mm)

Library info		Product info							
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C		
System	Kit								
AnyOne External Small Custom Abutment (9, 13mm)	ZrGEN 4015	Small	Ø4.2	0.6	4.5	Hex	AEXEPS 4015.MTN		
	ZrGEN 4016				6.0		AEXEPS 4016.MTN		
	ZrGEN 4018				8.0		AEXEPS 4018.MTN		
	ZrGEN 4025			1.5	4.5		AEXEPS 4025.MTN		
	ZrGEN 4026				6.0		AEXEPS 4026.MTN		
	ZrGEN 4028				8.0		AEXEPS 4028.MTN		
	ZrGEN 4035			3.0	4.5		AEXEPS 4035.MTN		
	ZrGEN 4036				6.0		AEXEPS 4036.MTN		
	ZrGEN 4038				8.0		AEXEPS 4038.MTN		
	ZrGEN 4045			4.0	4.5		AEXEPS 4045.MTN		
	ZrGEN 4046				6.0		AEXEPS 4046.MTN		
	ZrGEN 4048				8.0		AEXEPS 4048.MTN		
	ZrGEN 4515		Ø4.5	0.6	4.5	AEXEPS 4515.MTN			
	ZrGEN 4516				6.0	AEXEPS 4516.MTN			
	ZrGEN 4518				8.0	AEXEPS 4518.MTN			
	ZrGEN 4525			1.5	4.5	AEXEPS 4525.MTN			
	ZrGEN 4526				6.0	AEXEPS 4526.MTN			
	ZrGEN 4528				8.0	AEXEPS 4528.MTN			
	ZrGEN 4535			3.0	4.5	AEXEPS 4535.MTN			
	ZrGEN 4536				6.0	AEXEPS 4536.MTN			
	ZrGEN 4538				8.0	AEXEPS 4538.MTN			
	ZrGEN 4545			4.0	4.5	AEXEPS 4545.MTN			
	ZrGEN 4546				6.0	AEXEPS 4546.MTN			
	ZrGEN 4548				8.0	AEXEPS 4548.MTN			
AnyOne External Regular Custom Abutment (9, 13mm)	ZrGEN 4515	Regular	Ø4.5	0.6	4.5	Hex	AEXEPR 4515.MTN		
	ZrGEN 4516				6.0		AEXEPR 4516.MTN		
	ZrGEN 4518				8.0		AEXEPR 4518.MTN		
	ZrGEN 4525			1.5	4.5		AEXEPR 4525.MTN		
	ZrGEN 4526				6.0		AEXEPR 4526.MTN		
	ZrGEN 4528				8.0		AEXEPR 4528.MTN		
	ZrGEN 4535			3.0	4.5		AEXEPR 4535.MTN		
	ZrGEN 4536				6.0		AEXEPR 4536.MTN		
	ZrGEN 4538				8.0		AEXEPR 4538.MTN		
	ZrGEN 4545			4.0	4.5		AEXEPR 4545.MTN		
	ZrGEN 4546				6.0		AEXEPR 4546.MTN		
	ZrGEN 4548				8.0		AEXEPR 4548.MTN		
	AnyOne External Wide Custom Abutment (9, 13mm)		ZrGEN 5515	Wide	Ø5.5	0.6	4.5	Hex	AEXEPR 5515.MTN
			ZrGEN 5516				6.0		AEXEPR 5516.MTN
			ZrGEN 5518				8.0		AEXEPR 5518.MTN
			ZrGEN 5525			1.5	4.5		AEXEPR 5525.MTN
			ZrGEN 5526				6.0		AEXEPR 5526.MTN
			ZrGEN 5528				8.0		AEXEPR 5528.MTN
			ZrGEN 5535		3.0	4.5	AEXEPR 5535.MTN		
			ZrGEN 5536			6.0	AEXEPR 5536.MTN		
			ZrGEN 5538			8.0	AEXEPR 5538.MTN		
			ZrGEN 5545		4.0	4.5	AEXEPR 5545.MTN		
			ZrGEN 5546			6.0	AEXEPR 5546.MTN		
			ZrGEN 5548			8.0	AEXEPR 5548.MTN		



**TiGEN**  
Pre-milled



Library info		Product info				
Pre-milled Option In exocad library		Size	Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit					
 AnyOne External TiGEN Custom abutment	TiGEN 1220	Small	12	20	Hex	AETS1220.MTN
	TiGEN 1220	Regular	12	20		AETR1220.MTN
	TiGEN 1220	Wide	12	20		AETW1220.MTN



**ZrGEN**  
Ti-Base



AnyOne External  
Small Custom  
Abutment  
(9, 13mm)



AnyOne External  
Regular Custom  
Abutment  
(9, 13mm)

Library info			Product info					
ZrGEN Option In 3shape library			Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm						
	ZrGEN 4015_S09	ZrGEN 4015_S13	Small	Ø4.2	0.6	4.5	Hex	AEXEPS 4015.MTN
	ZrGEN 4016_S09	ZrGEN 4016_S13				6.0		AEXEPS 4016.MTN
	ZrGEN 4018_S09	ZrGEN 4018_S13				8.0		AEXEPS 4018.MTN
	ZrGEN 4025_S09	ZrGEN 4025_S13			1.5	4.5		AEXEPS 4025.MTN
	ZrGEN 4026_S09	ZrGEN 4026_S13				6.0		AEXEPS 4026.MTN
	ZrGEN 4028_S09	ZrGEN 4028_S13				8.0		AEXEPS 4028.MTN
	ZrGEN 4035_S09	ZrGEN 4035_S13			3.0	4.5		AEXEPS 4035.MTN
	ZrGEN 4036_S09	ZrGEN 4036_S13				6.0		AEXEPS 4036.MTN
	ZrGEN 4038_S09	ZrGEN 4038_S13				8.0		AEXEPS 4038.MTN
	ZrGEN 4045_S09	ZrGEN 4045_S13			4.0	4.5		AEXEPS 4045.MTN
	ZrGEN 4046_S09	ZrGEN 4046_S13				6.0		AEXEPS 4046.MTN
	ZrGEN 4048_S09	ZrGEN 4048_S13				8.0		AEXEPS 4048.MTN
	ZrGEN 4515_S09	ZrGEN 4515_S13		Ø4.5	0.6	4.5	AEXEPS 4515.MTN	
	ZrGEN 4516_S09	ZrGEN 4516_S13				6.0	AEXEPS 4516.MTN	
	ZrGEN 4518_S09	ZrGEN 4518_S13				8.0	AEXEPS 4518.MTN	
	ZrGEN 4525_S09	ZrGEN 4525_S13			1.5	4.5	AEXEPS 4525.MTN	
	ZrGEN 4526_S09	ZrGEN 4526_S13				6.0	AEXEPS 4526.MTN	
	ZrGEN 4528_S09	ZrGEN 4528_S13				8.0	AEXEPS 4528.MTN	
	ZrGEN 4535_S09	ZrGEN 4535_S13			3.0	4.5	AEXEPS 4535.MTN	
	ZrGEN 4536_S09	ZrGEN 4536_S13				6.0	AEXEPS 4536.MTN	
	ZrGEN 4538_S09	ZrGEN 4538_S13				8.0	AEXEPS 4538.MTN	
	ZrGEN 4545_S09	ZrGEN 4545_S13			4.0	4.5	AEXEPS 4545.MTN	
	ZrGEN 4546_S09	ZrGEN 4546_S13				6.0	AEXEPS 4546.MTN	
	ZrGEN 4548_S09	ZrGEN 4548_S13				8.0	AEXEPS 4548.MTN	
	ZrGEN 4515_S09	ZrGEN 4515_S13	Regular	Ø4.5	0.6	4.5	Hex	AEXEPR 4515.MTN
	ZrGEN 4516_S09	ZrGEN 4516_S13				6.0		AEXEPR 4516.MTN
	ZrGEN 4518_S09	ZrGEN 4518_S13				8.0		AEXEPR 4518.MTN
	ZrGEN 4525_S09	ZrGEN 4525_S13			1.5	4.5		AEXEPR 4525.MTN
	ZrGEN 4526_S09	ZrGEN 4526_S13				6.0		AEXEPR 4526.MTN
	ZrGEN 4528_S09	ZrGEN 4528_S13				8.0		AEXEPR 4528.MTN
	ZrGEN 4535_S09	ZrGEN 4535_S13			3.0	4.5		AEXEPR 4535.MTN
	ZrGEN 4536_S09	ZrGEN 4536_S13				6.0		AEXEPR 4536.MTN
	ZrGEN 4538_S09	ZrGEN 4538_S13				8.0		AEXEPR 4538.MTN
	ZrGEN 4545_S09	ZrGEN 4545_S13			4.0	4.5		AEXEPR 4545.MTN
	ZrGEN 4546_S09	ZrGEN 4546_S13				6.0		AEXEPR 4546.MTN
	ZrGEN 4548_S09	ZrGEN 4548_S13				8.0		AEXEPR 4548.MTN
	ZrGEN 5515_S09	ZrGEN 5515_S13	Wide	Ø5.5	0.6	4.5	Hex	AEXEPR 5515.MTN
	ZrGEN 5516_S09	ZrGEN 5516_S13				6.0		AEXEPR 5516.MTN
	ZrGEN 5518_S09	ZrGEN 5518_S13				8.0		AEXEPR 5518.MTN
	ZrGEN 5525_S09	ZrGEN 5525_S13			1.5	4.5		AEXEPR 5525.MTN
	ZrGEN 5526_S09	ZrGEN 5526_S13				6.0		AEXEPR 5526.MTN
	ZrGEN 5528_S09	ZrGEN 5528_S13				8.0		AEXEPR 5528.MTN
	ZrGEN 5535_S09	ZrGEN 5535_S13			3.0	4.5		AEXEPR 5535.MTN
	ZrGEN 5536_S09	ZrGEN 5536_S13				6.0		AEXEPR 5536.MTN
	ZrGEN 5538_S09	ZrGEN 5538_S13				8.0		AEXEPR 5538.MTN
	ZrGEN 5545_S09	ZrGEN 5545_S13			4.0	4.5		AEXEPR 5545.MTN
	ZrGEN 5546_S09	ZrGEN 5546_S13				6.0		AEXEPR 5546.MTN
	ZrGEN 5548_S09	ZrGEN 5548_S13				8.0		AEXEPR 5548.MTN






**TiGEN**  
Pre-milled



AnyOne\_External\_TiGEN  
Custom abutment


Library info			Product info				
Pre-milled Option In 3shape library			Size	Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
AnyOne_External_TiGEN Custom abutment	TiGEN 1220_S09	TiGEN 1220_S13	Small	12	20	Hex	AETS1220.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	Regular	12	20		AETRI220.MTN
	TiGEN 1220_S09	TiGEN 1220_S13	Wide	12	20		AETW1220.MTN





\* EZ Post Prosthesis by utilizing the stock abutment like a scan abutment.

# Scan abutment



Scan Abutment [MN]

Product info			Library info (Choose below CAD library option)	
Diameter	Lenght	Ref.C	exocad library name	3shape library name
Ø4.0	9	MISS3509T	MiNi_Custom Abutment (9mm)	MiNi_ZrGEN_Custom abutment
	13	MISS3513T	MiNi_Custom Abutment (13mm)	MINI_TiGEN_Custom abutment



**ZrGEN**

Ti-Base



Mini Custom Abutment  
(9, 13mm)

Library info		Product info				
ZrGEN Option In exocad library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
	ZrGEN 3013	Ø3.0	0.6	2.5	Hex	MIPN3013.MTN
	ZrGEN 3013	Ø3.0	0.6	2.5	Non-Hex	MIPN3013N.MTN



**TiGEN**  
Pre-milled





Mini Custom Abutment  
(9mm, 13mm)

Library info		Product info			
Pre-milled Option In exocad library		Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit				
	TiGEN 1020	10	20	Hex	MITN1020.MTN
	TiGEN 1020	10	20	Non-Hex	MITN1020N.MTN



**EZ Post**  
Stock abutment




Library info		Product info				
EZ Post Option In exocad library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
 MiNi Scanbody - EZpost (9, 13mm)	EZ Post 3505	Ø3.5	1.0	5.5	Hex	MIEP3505HT
	EZ Post 3515		1.5			MIEP3515HT
	EZ Post 3525		2.5			MIEP3525HT
	EZ Post 3535		3.5			MIEP3535HT
	EZ Post 3545		4.5			MIEP3545HT
	EZ Post 3507		1.0	7.0		MIEP3507HT
	EZ Post 3517		1.5			MIEP3517HT
	EZ Post 3527		2.5			MIEP3527HT
	EZ Post 3537		3.5			MIEP3537HT
	EZ Post 3547		4.5			MIEP3547HT
 MiNi EZpost - EZpost	EZ Post 3509	Ø3.5	1.0	9.0	Hex	MIEP3509HT
	EZ Post 3519		1.5			MIEP3519HT
	EZ Post 3529		2.5			MIEP3529HT
	EZ Post 3539		3.5			MIEP3539HT
	EZ Post 3549		4.5			MIEP3549HT



**ZrGEN**

Ti-Base



Library info			Product info				
ZrGEN Option In 3shape library			Diameter (∅)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
	ZrGEN 3013_S09	ZrGEN 3013_S13	∅3.0	0.6	4.5	Hex	MIPN3013.MTN
 Mini ZrGEN Custom Abutment (9, 13mm)	ZrGEN 3013_S09	ZrGEN 3013.S13	∅3.0	0.6	4.5	Non- Hex	AEXEPR 5548.MTN



**TiGEN**  
Pre-milled



Mini TiGEN  
Custom abutment  
(9, 13mm)


Library info			Product info				
Pre-milled Option In 3shape library			Size	Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
Mini TiGEN Custom abutment (9, 13mm)	TiGEN 1020_S09	TiGEN 1020_S13	Small	12	20	Hex	MITN1020.MTN
	TiGEN 1020_S09	TiGEN 1020_S13	Small	12	20	Non-hex	MITN1020N.MTN






**EZ Post**  
Stock abutment



Library info		Product info				
EZ Post Option In 3shape library		Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit					
	EZ Post 3505	Ø3.5	1.0	5.5	Hex	MIEP3505HT
	EZ Post 3515		1.5			MIEP3515HT
	EZ Post 3525		2.5			MIEP3525HT
	EZ Post 3535		3.5			MIEP3535HT
	EZ Post 3545		4.5			MIEP3545HT
 MiNi Scanbody - EZpost (9, 11mm)	EZ Post 3507		1.0	7.0		MIEP3507HT
	EZ Post 3517		1.5			MIEP3517HT
	EZ Post 3527		2.5			MIEP3527HT
	EZ Post 3537		3.5			MIEP3537HT
	EZ Post 3547		4.5			MIEP3547HT
 MiNi EZpost - EZpost	EZ Post 3509		1.0	9.0		MIEP3509HT
	EZ Post 3519		1.5			MIEP3519HT
	EZ Post 3529		2.5			MIEP3529HT
	EZ Post 3539		3.5			MIEP3539HT
	EZ Post 3549		4.5			MIEP3549HT




ST



\* EZ Post Prosthesis by utilizing the stock abutment like a scan abutment.

# Scan abutment

## Scan Abutment [ST]



Product info				Library info (Choose below CAD library option)	
Diameter	Lenght	Type	Ref.C	exocad library name	3shape library name
Ø4.0	9	Small	OSGSSC3110T	ST_Small Custom Abutment (9mm)	ST_Small_ZrGEN_Custom abutment ST_Regular_ZrGEN_Custom abutment ST_TiGEN_Custom abutment
	13		OSGSSC3111T	ST_Small Custom Abutment (13mm)	
	9	Regular	OSGSSC3210T	ST_Regular Custom Abutment (9mm)	
	13		OSGSSC3211T	ST_Regular Custom Abutment (13mm)	



**ZrGEN**

Ti-Base



ST\_Small Custom Abutment (9, 13mm)

Library info		Product info								
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C			
System	Kit									
ST_Small Custom Abutment (9, 13mm)	ZrGEN 3111	Small	Ø4.0	0.6	4.5	Hex	OSGSPA3111.MTN			
	ZrGEN 3112				6.0		OSGSPA3112.MTN			
	ZrGEN 3113				8.0		OSGSPA3113.MTN			
	ZrGEN 3121			1.5	4.5		OSGSPA3121.MTN			
	ZrGEN 3122				6.0		OSGSPA3122.MTN			
	ZrGEN 3123				8.0		OSGSPA3231.MTN			
	ZrGEN 3131			3.0	4.5		OSGSPA3131.MTN			
	ZrGEN 3132				6.0		OSGSPA3132.MTN			
	ZrGEN 3133				8.0		OSGSPA3133.MTN			
	ZrGEN 3141			4.0	4.5		OSGSPA3141.MTN			
	ZrGEN 3142				6.0		OSGSPA3142.MTN			
	ZrGEN 3143				8.0		OSGSPA3143.MTN			
	ZrGEN 3111			Small	Ø4.0		0.6	4.5	Non-Hex	OSGSPA3111N.MTN
	ZrGEN 3112							6.0		OSGSPA3112N.MTN
ZrGEN 3113	8.0	OSGSPA3113N.MTN								
ZrGEN 3121	1.5	4.5	OSGSPA3121N.MTN							
ZrGEN 3122		6.0	OSGSPA3122N.MTN							
ZrGEN 3123		8.0	OSGSPA3123N.MTN							
ZrGEN 3131	3.0	4.5	OSGSPA3131N.MTN							
ZrGEN 3132		6.0	OSGSPA3132N.MTN							
ZrGEN 3133		8.0	OSGSPA3133N.MTN							
ZrGEN 3141	4.0	4.5	OSGSPA3141N.MTN							
ZrGEN 3142		6.0	OSGSPA3142N.MTN							
ZrGEN 3143		8.0	OSGSPA3143N.MTN							

**ZrGEN**  
Ti-Base



Library info		Product info					
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit						
ST_Regular Custom Abutment (9, 13mm)	ZrGEN 3211	Regular	Ø4.0	0.6	4.5	Hex	OSGSPA3211.MTN
	ZrGEN 3212				6.0		OSGSPA3212.MTN
	ZrGEN 4018				8.0		OSGSPA4018.MTN
	ZrGEN 4025			1.5	4.5		OSGSPA4025.MTN
	ZrGEN 4026				6.0		OSGSPA4026.MTN
	ZrGEN 4028				8.0		OSGSPA4028.MTN
	ZrGEN 4035			3.0	4.5		OSGSPA4035.MTN
	ZrGEN 4036				6.0		OSGSPA4036.MTN
	ZrGEN 4038				8.0		OSGSPA4038.MTN
	ZrGEN 4045			4.0	4.5		OSGSPA4045.MTN
	ZrGEN 4046				6.0		OSGSPA4046.MTN
	ZrGEN 4048				8.0		OSGSPA4048.MTN
	ZrGEN 3211		Ø4.0	0.6	4.5	Non-Hex	OSGSPA3211N.MTN
	ZrGEN 3212				6.0		OSGSPA3212N.MTN
	ZrGEN 4018				8.0		OSGSPA4018N.MTN
	ZrGEN 4025			1.5	4.5		OSGSPA4025N.MTN
	ZrGEN 4026				6.0		OSGSPA4026N.MTN
	ZrGEN 4028				8.0		OSGSPA4028N.MTN
	ZrGEN 4035			3.0	4.5		OSGSPA4035N.MTN
	ZrGEN 4036				6.0		OSGSPA4036N.MTN
	ZrGEN 4038				8.0		OSGSPA4038N.MTN
	ZrGEN 4045			4.0	4.5		OSGSPA4045N.MTN
	ZrGEN 4046				6.0		OSGSPA4046N.MTN
	ZrGEN 4048				8.0		OSGSPA4048N.MTN
ZrGEN 4515	Ø4.5	0.6	4.5	Hex	OSGSPA4515.MTN		
ZrGEN 4516			6.0		OSGSPA4516.MTN		
ZrGEN 4518			8.0		OSGSPA4518.MTN		
ZrGEN 3221		1.5	4.5		OSGSPA3221.MTN		
ZrGEN 3222			6.0		OSGSPA3222.MTN		
ZrGEN 4528			8.0		OSGSPA4528.MTN		
ZrGEN 4535		3.0	4.5		OSGSPA4535.MTN		
ZrGEN 4536			6.0		OSGSPA4536.MTN		
ZrGEN 4538			8.0		OSGSPA4538.MTN		
ZrGEN 4545		4.0	4.5		OSGSPA4545.MTN		
ZrGEN 4546			6.0		OSGSPA4546.MTN		
ZrGEN 4548			8.0		OSGSPA4548.MTN		
ZrGEN 4515	Ø4.5	0.6	4.5	Non-Hex	OSGSPA4515N.MTN		
ZrGEN 4516			6.0		OSGSPA4516N.MTN		
ZrGEN 4518			8.0		OSGSPA4518N.MTN		
ZrGEN 3221		1.5	4.5		OSGSPA3221N.MTN		
ZrGEN 3222			6.0		OSGSPA3222N.MTN		
ZrGEN 4528			8.0		OSGSPA4528N.MTN		
ZrGEN 4535		3.0	4.5		OSGSPA4535N.MTN		
ZrGEN 4536			6.0		OSGSPA4536N.MTN		
ZrGEN 4538			8.0		OSGSPA4538N.MTN		
ZrGEN 4545		4.0	4.5		OSGSPA4545N.MTN		
ZrGEN 4546			6.0		OSGSPA4546N.MTN		
ZrGEN 4548			8.0		OSGSPA4548N.MTN		



**TiGEN**  
Pre-milled



ST\_Small Custom Abutment  
ST\_Regular Custom Abutment  
(9, 13mm)


Library info		Product info				
Pre-milled Option In exocad library		Size	Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit					
	TiGEN 1020	Small	10	20	Hex	OSTG3112.MTN
	TiGEN 1220		12	20		OSTG3111.MTN
	TiGEN 1020	Regular	10	20		OSTG3212.MTN
	TiGEN 1220		12	20		OSTG3211.MTN
	TiGEN 1020	Small	10	20	Non-Hex	OSTG3112N.MTN
	TiGEN 1220		12	20		OSTG3111N.MTN
	TiGEN 1020	Regular	10	20		OSTG3212N.MTN
	TiGEN 1220		12	20		OSTG3211N.MTN



**EZ Post**  
Stock abutment




Library info		Product info								
EZ Post Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C			
System	Kit									
ST Scanbody - EZpost (9, 13mm)	EZ Post S461I	Small	Ø4.5	1.0	5.5	Hex	STEA461IT			
	EZ Post S462I			2.0			STEA462IT			
	EZ Post S463I			3.0			STEA463IT			
	EZ Post S464I			4.0			STEA464IT			
	EZ Post S465I			5.0			STEA465IT			
	EZ Post S471I			1.0	7.0		STEA471IT			
	EZ Post S472I			2.0			STEA472IT			
	EZ Post S473I			3.0			STEA473IT			
	EZ Post S474I			4.0			STEA474IT			
	EZ Post S475I			5.0			STEA475IT			
	EZ Post S461I			1.0	Non-Hex		STEA461INT			
	EZ Post S462I			2.0			STEA462INT			
	EZ Post S463I			3.0			STEA463INT			
	EZ Post S464I			4.0			STEA464INT			
	EZ Post S465I			5.0			STEA465INT			
	EZ Post S471I		1.0	7.0		STEA471INT				
	EZ Post S472I		2.0			STEA472INT				
	EZ Post S473I		3.0			STEA473INT				
	EZ Post S474I		4.0			STEA474INT				
	EZ Post S475I		5.0			STEA475INT				
	ST EZpost - EZpost		EZ Post 461I	Regular		Ø4.5	1.0	5.5	Hex	STEAS461IT
			EZ Post 462I				2.0			STEAS462IT
			EZ Post 463I				3.0			STEAS463IT
			EZ Post 464I				4.0			STEAS464IT
			EZ Post 465I				5.0			STEAS465IT
			EZ Post 471I		1.0		7.0	STEAS471IT		
			EZ Post 472I		2.0			STEAS472IT		
			EZ Post 473I		3.0			STEAS473IT		
			EZ Post 474I		4.0			STEAS474IT		
			EZ Post 475I		5.0			STEAS475IT		
EZ Post 461I		1.0	Non-Hex		STEAS461INT					
EZ Post 462I		2.0			STEAS462INT					
EZ Post 463I		3.0			STEAS463INT					
EZ Post 464I		4.0			STEAS464INT					
EZ Post 465I		5.0			STEAS465INT					
EZ Post 471I		1.0			7.0	STEAS471INT				
EZ Post 472I		2.0				STEAS472INT				
EZ Post 473I		3.0				STEAS473INT				
EZ Post 474I		4.0				STEAS474INT				
EZ Post 475I		5.0				STEAS475INT				
EZ Post 5410		1.0			Ø5.0	4.0	Hex	STEA5410T		
EZ Post 5420		2.0						STEA5420T		
EZ Post 5430		3.0						STEA5430T		
EZ Post 5440		4.0						STEA5440T		
EZ Post 5450		5.0						STEA5450T		
EZ Post 5610		1.0	5.5			Hex		STEA5610T		
EZ Post 5620		2.0						STEA5620T		
EZ Post 5630		3.0						STEA5630T		
EZ Post 5640		4.0						STEA5640T		
EZ Post 5650		5.0						STEA5650T		
EZ Post 5710	1.0	7.0	Hex	STEA5710T						
EZ Post 5720	2.0			STEA5720T						
EZ Post 5730	3.0			STEA5730T						
EZ Post 5740	4.0			STEA5740T						
EZ Post 5750	5.0			STEA5750T						
EZ Post 5410	1.0	Non-Hex	4.0	Non-Hex	STEA5410NT					
EZ Post 5420	2.0				STEA5420NT					
EZ Post 5430	3.0				STEA5430NT					
EZ Post 5440	4.0				STEA5440NT					
EZ Post 5450	5.0				STEA5450NT					
EZ Post 5610	1.0		5.5		Non-Hex	STEA5610NT				
EZ Post 5620	2.0					STEA5620NT				
EZ Post 5630	3.0					STEA5630NT				
EZ Post 5640	4.0					STEA5640NT				
EZ Post 5650	5.0					STEA5650NT				
EZ Post 5710	1.0		7.0		Non-Hex	STEA5710NT				
EZ Post 5720	2.0					STEA5720NT				
EZ Post 5730	3.0					STEA5730NT				
EZ Post 5740	4.0					STEA5740NT				
EZ Post 5750	5.0					STEA5750NT				



**EZ Post**  
Stock abutment



Library info		Product info					
EZ Post Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit						
	EZ Post 6410	Regular	Ø6.0	1.0	4.0	Hex	STEA6410T
	EZ Post 6420			2.0			STEA6420T
	EZ Post 6430			3.0			STEA6430T
	EZ Post 6440			4.0			STEA6440T
	EZ Post 6450			5.0			STEA6450T
	EZ Post 6610			1.0	5.5		STEA6610T
	EZ Post 6620			2.0			STEA6620T
	EZ Post 6630			3.0			STEA6630T
	EZ Post 6640			4.0			STEA6640T
	EZ Post 6650			5.0			STEA6650T
	EZ Post 6710			1.0	7.0		STEA6710T
	EZ Post 6720			2.0			STEA6720T
	EZ Post 6730			3.0			STEA6730T
	EZ Post 6740			4.0			STEA6740T
	EZ Post 6750			5.0			STEA6750T
	EZ Post 6410		Regular	1.0	4.0	Hex	STEA6410NT
	EZ Post 6420	2.0		STEA6420NT			
	EZ Post 6430	3.0		STEA6430NT			
	EZ Post 6440	4.0		STEA6440NT			
	EZ Post 6450	5.0		STEA6450NT			
	EZ Post 6610	1.0		5.5	STEA6610NT		
	EZ Post 6620	2.0			STEA6620NT		
	EZ Post 6630	3.0			STEA6630NT		
	EZ Post 6640	4.0			STEA6640NT		
	EZ Post 6650	5.0			STEA6650NT		
	EZ Post 6710	1.0		7.0	STEA6710NT		
	EZ Post 6720	2.0			STEA6720NT		
	EZ Post 6730	3.0			STEA6730NT		
	EZ Post 6740	4.0			STEA6740NT		
	EZ Post 6750	5.0			STEA6750NT		
	EZ Post 7610		Ø7.0	1.0	5.5	Hex	STEA7610T
	EZ Post 7620	2.0		STEA7620T			
	EZ Post 7630	3.0		STEA7630T			
	EZ Post 7640	4.0		STEA7640T			
	EZ Post 7650	5.0		STEA7650T			
	EZ Post 7610		Ø7.0	1.0	5.5	Non-Hex	STEA7610NT
	EZ Post 7620	2.0		STEA7620NT			
	EZ Post 7630	3.0		STEA7630NT			
	EZ Post 7640	4.0		STEA7640NT			
	EZ Post 7650	5.0		STEA7650NT			





**ZrGEN**  
Ti-Base



Library info			Product info					
ZrGEN Option In 3shape library			Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm						
ST_ZrGEN Custom abutment	ZrGEN 3111_S09	ZrGEN 3111_S13	Small	Ø4.0	0.6	4.5	Hex	OSGSPA3111.MTN
	ZrGEN 3112_S09	ZrGEN 3112_S13				6.0		OSGSPA3112.MTN
	ZrGEN 3113_S09	ZrGEN 3113_S13				8.0		OSGSPA3113.MTN
	ZrGEN 3121_S09	ZrGEN 3121_S13			1.5	4.5		OSGSPA3121.MTN
	ZrGEN 3122_S09	ZrGEN 3122_S13				6.0		OSGSPA3122.MTN
	ZrGEN 3123_S09	ZrGEN 3123_S13				8.0		OSGSPA3123.MTN
	ZrGEN 3131_S09	ZrGEN 3131_S13			3.0	4.5		OSGSPA3131.MTN
	ZrGEN 3132_S09	ZrGEN 3132_S13				6.0		OSGSPA3132.MTN
	ZrGEN 3133_S09	ZrGEN 3133_S13				8.0		OSGSPA3133.MTN
	ZrGEN 3141_S09	ZrGEN 3141_S13			4.0	4.5		OSGSPA3141.MTN
	ZrGEN 3142_S09	ZrGEN 3142_S13				6.0		OSGSPA3142.MTN
	ZrGEN 3143_S09	ZrGEN 3143_S13				8.0		OSGSPA3143.MTN
	ZrGEN 3111_S09	ZrGEN 3111_S13			Small	Ø4.0		0.6
ZrGEN 3112_S09	ZrGEN 3112_S13	6.0	OSGSPA3112N.MTN					
ZrGEN 3113_S09	ZrGEN 3113_S13	8.0	OSGSPA3113N.MTN					
ZrGEN 3121_S09	ZrGEN 3121_S13	1.5	4.5	OSGSPA3121N.MTN				
ZrGEN 3122_S09	ZrGEN 3122_S13		6.0	OSGSPA3122N.MTN				
ZrGEN 3123_S09	ZrGEN 3123_S13		8.0	OSGSPA3123N.MTN				
ZrGEN 3131_S09	ZrGEN 3131_S13	3.0	4.5	OSGSPA3131N.MTN				
ZrGEN 3132_S09	ZrGEN 3132_S13		6.0	OSGSPA3132N.MTN				
ZrGEN 3133_S09	ZrGEN 3133_S13		8.0	OSGSPA3133N.MTN				
ZrGEN 3141_S09	ZrGEN 3141_S13	4.0	4.5	OSGSPA3141N.MTN				
ZrGEN 3142_S09	ZrGEN 3142_S13		6.0	OSGSPA3142N.MTN				
ZrGEN 3143_S09	ZrGEN 3143_S13		8.0	OSGSPA3143N.MTN				



**ZrGEN**

Ti-Base



ST\_ZrGEN  
Custom Abutment


Library info			Product info					
ZrGEN Option In 3shape library			Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm						
	ZrGEN 3211_S09	ZrGEN 3211_S13	Regular	Ø4.0	0.6	4.5	Hex	OSGSPA3211.MTN
	ZrGEN 3212_S09	ZrGEN 3212_S13				6.0		OSGSPA3212.MTN
	ZrGEN 4018_S09	ZrGEN 4018_S13				8.0		OSGSPA4018.MTN
	ZrGEN 4025_S09	ZrGEN 4025_S13			1.5	4.5		OSGSPA4025.MTN
	ZrGEN 4026_S09	ZrGEN 4026_S13				6.0		OSGSPA4026.MTN
	ZrGEN 4028_S09	ZrGEN 4028_S13				8.0		OSGSPA4028.MTN
	ZrGEN 4035_S09	ZrGEN 4035_S13			3.0	4.5		OSGSPA4035.MTN
	ZrGEN 4036_S09	ZrGEN 4036_S13				6.0		OSGSPA4036.MTN
	ZrGEN 4038_S09	ZrGEN 4038_S13				8.0		OSGSPA4038.MTN
	ZrGEN 4045_S09	ZrGEN 4045_S13			4.0	4.5		OSGSPA4045.MTN
	ZrGEN 4046_S09	ZrGEN 4046_S13				6.0		OSGSPA4046.MTN
	ZrGEN 4048_S09	ZrGEN 4048_S13				8.0		OSGSPA4048.MTN
	ZrGEN 3211_S09	ZrGEN 3211_S13		Ø4.0	0.6	4.5	Non-Hex	OSGSPA3211N.MTN
	ZrGEN 3212_S09	ZrGEN 3212_S13				6.0		OSGSPA3212N.MTN
	ZrGEN 4018_S09	ZrGEN 4018_S13				8.0		OSGSPA4018N.MTN
	ZrGEN 4025_S09	ZrGEN 4025_S13			1.5	4.5		OSGSPA4025N.MTN
	ZrGEN 4026_S09	ZrGEN 4026_S13				6.0		OSGSPA4026N.MTN
	ZrGEN 4028_S09	ZrGEN 4028_S13				8.0		OSGSPA4028N.MTN
	ZrGEN 4035_S09	ZrGEN 4035_S13			3.0	4.5		OSGSPA4035N.MTN
	ZrGEN 4036_S09	ZrGEN 4036_S13				6.0		OSGSPA4036N.MTN
	ZrGEN 4038_S09	ZrGEN 4038_S13				8.0		OSGSPA4038N.MTN
	ZrGEN 4045_S09	ZrGEN 4045_S13			4.0	4.5		OSGSPA4045N.MTN
	ZrGEN 4046_S09	ZrGEN 4046_S13				6.0		OSGSPA4046N.MTN
	ZrGEN 4048_S09	ZrGEN 4048_S13				8.0		OSGSPA4048N.MTN
	ZrGEN 4515_S09	ZrGEN 4515_S13	Ø4.5	0.6	4.5	Hex	OSGSPA4515.MTN	
	ZrGEN 4516_S09	ZrGEN 4516_S13			6.0		OSGSPA4516.MTN	
	ZrGEN 4518_S09	ZrGEN 4518_S13			8.0		OSGSPA4518.MTN	
	ZrGEN 3221_S09	ZrGEN 3221_S13		1.5	4.5		OSGSPA3221.MTN	
	ZrGEN 3222_S09	ZrGEN 3222_S13			6.0		OSGSPA3222.MTN	
	ZrGEN 4528_S09	ZrGEN 4528_S13			8.0		OSGSPA4528.MTN	
	ZrGEN 4535_S09	ZrGEN 4535_S13		3.0	4.5		OSGSPA4535.MTN	
	ZrGEN 4536_S09	ZrGEN 4536_S13			6.0		OSGSPA4536.MTN	
	ZrGEN 4538_S09	ZrGEN 4538_S13			8.0		OSGSPA4538.MTN	
	ZrGEN 4545_S09	ZrGEN 4545_S13		4.0	4.5		OSGSPA4545.MTN	
	ZrGEN 4546_S09	ZrGEN 4546_S13			6.0		OSGSPA4546.MTN	
	ZrGEN 4548_S09	ZrGEN 4548_S13			8.0		OSGSPA4548.MTN	
	ZrGEN 4515_S09	ZrGEN 4515_S13	Ø4.5	0.6	4.5	Non-Hex	OSGSPA4515N.MTN	
	ZrGEN 4516_S09	ZrGEN 4516_S13			6.0		OSGSPA4516N.MTN	
	ZrGEN 4518_S09	ZrGEN 4518_S13			8.0		OSGSPA4518N.MTN	
	ZrGEN 3221_S09	ZrGEN 3221_S13		1.5	4.5		OSGSPA3221N.MTN	
	ZrGEN 3222_S09	ZrGEN 3222_S13			6.0		OSGSPA3222N.MTN	
	ZrGEN 4528_S09	ZrGEN 4528_S13			8.0		OSGSPA4528N.MTN	
	ZrGEN 4535_S09	ZrGEN 4535_S13		3.0	4.5		OSGSPA4535N.MTN	
	ZrGEN 4536_S09	ZrGEN 4536_S13			6.0		OSGSPA4536N.MTN	
	ZrGEN 4538_S09	ZrGEN 4538_S13			8.0		OSGSPA4538N.MTN	
	ZrGEN 4545_S09	ZrGEN 4545_S13		4.0	4.5		OSGSPA4545N.MTN	
	ZrGEN 4546_S09	ZrGEN 4546_S13			6.0		OSGSPA4546N.MTN	
	ZrGEN 4548_S09	ZrGEN 4548_S13			8.0		OSGSPA4548N.MTN	



**TiGEN**  
Pre-milled



ST\_TiGEN  
Custom abutment


Library info			Product info				
Pre-milled Option In 3shape library			Size	Diameter (∅)	Length (mm)	Type	Ref.C
System	Kit for 9mm	Kit for 13mm					
 ST_TiGEN Custom abutment	TiGEN 1020_S09	TiGEN 1020_S13	Small	10	20	Hex	OSTG3112.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		OSTG3111.MTN
	TiGEN 1020_S09	TiGEN 1020_S13	Regular	10	20		OSTG3212.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		OSTG3211.MTN
	TiGEN 1020_S09	TiGEN 1020_S13	Small	10	20	Non-hex	OSTG3112N.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		OSTG3111N.MTN
	TiGEN 1020_S09	TiGEN 1020_S13	Regular	10	20		OSTG3212N.MTN
	TiGEN 1220_S09	TiGEN 1220_S13		12	20		OSTG3211N.MTN



**EZ Post**  
Stock abutment



Library info		Product info							
EZ Post Option In 3shape library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C		
System	Kit								
ST_EZ Post_Scan Abutment	EZ Post S4611	Small	Ø4.5	1.0	5.5	Hex	STEA4611T		
	EZ Post S4621			2.0			STEA4621T		
	EZ Post S4631			3.0			STEA4631T		
	EZ Post S4641			4.0			STEA4641T		
	EZ Post S4651			5.0			STEA4651T		
	EZ Post S4711			1.0	7.0		STEA4711T		
	EZ Post S4721			2.0			STEA4721T		
	EZ Post S4731			3.0			STEA4731T		
	EZ Post S4741			4.0			STEA4741T		
	EZ Post S4751			5.0			STEA4751T		
	EZ Post S4611		1.0	Non-Hex	STEA4611NT				
	EZ Post S4621		2.0		STEA4621NT				
	EZ Post S4631		3.0		STEA4631NT				
	EZ Post S4641		4.0		STEA4641NT				
	EZ Post S4651		5.0		STEA4651NT				
	EZ Post S4711		1.0		7.0	STEA4711NT			
	EZ Post S4721		2.0			STEA4721NT			
	EZ Post S4731		3.0			STEA4731NT			
	EZ Post S4741		4.0			STEA4741NT			
	EZ Post S4751		5.0			STEA4751NT			
	ST_EZ Post_EZ Post		EZ Post 4611	Regular	Ø4.5	1.0	5.5	Hex	STEAS4611T
			EZ Post 4621			2.0			STEAS4621T
			EZ Post 4631			3.0			STEAS4631T
			EZ Post 4641			4.0			STEAS4641T
			EZ Post 4651			5.0			STEAS4651T
			EZ Post 4711			1.0	7.0		STEAS4711T
			EZ Post 4721			2.0			STEAS4721T
			EZ Post 4731			3.0			STEAS4731T
			EZ Post 4741			4.0			STEAS4741T
			EZ Post 4751			5.0			STEAS4751T
			EZ Post 4611		1.0	Non-Hex	STEAS4611NT		
			EZ Post 4621		2.0		STEAS4621NT		
			EZ Post 4631		3.0		STEAS4631NT		
			EZ Post 4641		4.0		STEAS4641NT		
			EZ Post 4651		5.0		STEAS4651NT		
EZ Post 4711		1.0	7.0		STEAS4711NT				
EZ Post 4721		2.0			STEAS4721NT				
EZ Post 4731		3.0			STEAS4731NT				
EZ Post 4741		4.0			STEAS4741NT				
EZ Post 4751		5.0			STEAS4751NT				
EZ Post 5410		1.0	Ø5.0		4.0	Hex	STEA5410T		
EZ Post 5420		2.0					STEA5420T		
EZ Post 5430		3.0					STEA5430T		
EZ Post 5440		4.0					STEA5440T		
EZ Post 5450		5.0					STEA5450T		
EZ Post 5610		1.0			5.5		Hex	STEA5610T	
EZ Post 5620		2.0						STEA5620T	
EZ Post 5630		3.0						STEA5630T	
EZ Post 5640		4.0						STEA5640T	
EZ Post 5650		5.0						STEA5650T	
EZ Post 5710		1.0	7.0		Hex	STEA5710T			
EZ Post 5720		2.0				STEA5720T			
EZ Post 5730		3.0				STEA5730T			
EZ Post 5740		4.0				STEA5740T			
EZ Post 5750		5.0				STEA5750T			
EZ Post 5410	1.0	Non-Hex	4.0	Non-Hex		STEA5410NT			
EZ Post 5420	2.0					STEA5420NT			
EZ Post 5430	3.0					STEA5430NT			
EZ Post 5440	4.0					STEA5440NT			
EZ Post 5450	5.0					STEA5450NT			
EZ Post 5610	1.0		5.5		Non-Hex	STEA5610NT			
EZ Post 5620	2.0					STEA5620NT			
EZ Post 5630	3.0					STEA5630NT			
EZ Post 5640	4.0					STEA5640NT			
EZ Post 5650	5.0					STEA5650NT			
EZ Post 5710	1.0	7.0	Non-Hex	STEA5710NT					
EZ Post 5720	2.0			STEA5720NT					
EZ Post 5730	3.0			STEA5730NT					
EZ Post 5740	4.0			STEA5740NT					
EZ Post 5750	5.0			STEA5750NT					



**EZ Post**  
Stock abutment




Library info  
EZ Post Option In 3shape library

System Kit


Product info

Size Diameter (Ø) Cuff Height (mm) Post Height (mm) Type Ref.C



Library info		Product info							
EZ Post Option In 3shape library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C		
System	Kit								
	EZ Post 6410	Regular	Ø6.0	1.0	4.0	Hex	STEA6410T		
	EZ Post 6420			2.0			STEA6420T		
	EZ Post 6430			3.0			STEA6430T		
	EZ Post 6440			4.0			STEA6440T		
	EZ Post 6450			5.0			STEA6450T		
	EZ Post 6610			1.0	5.5		Hex	STEA6610T	
	EZ Post 6620			2.0				STEA6620T	
	EZ Post 6630			3.0				STEA6630T	
	EZ Post 6640			4.0				STEA6640T	
	EZ Post 6650			5.0				STEA6650T	
	EZ Post 6710			1.0	7.0			Hex	STEA6710T
	EZ Post 6720			2.0					STEA6720T
	EZ Post 6730			3.0					STEA6730T
	EZ Post 6740			4.0					STEA6740T
	EZ Post 6750			5.0					STEA6750T
	EZ Post 6410	1.0	4.0	Hex	STEA6410NT				
	EZ Post 6420	2.0			STEA6420NT				
	EZ Post 6430	3.0			STEA6430NT				
	EZ Post 6440	4.0			STEA6440NT				
	EZ Post 6450	5.0			STEA6450NT				
	EZ Post 6610	1.0	5.5		Hex	STEA6610NT			
	EZ Post 6620	2.0				STEA6620NT			
	EZ Post 6630	3.0				STEA6630NT			
	EZ Post 6640	4.0				STEA6640NT			
	EZ Post 6650	5.0				STEA6650NT			
	EZ Post 6710	1.0	7.0			Hex	STEA6710NT		
	EZ Post 6720	2.0					STEA6720NT		
	EZ Post 6730	3.0					STEA6730NT		
	EZ Post 6740	4.0					STEA6740NT		
	EZ Post 6750	5.0					STEA6750NT		
	EZ Post 7610	1.0	Ø7.0	5.5			Hex	STEA7610T	
	EZ Post 7620	2.0						STEA7620T	
	EZ Post 7630	3.0						STEA7630T	
	EZ Post 7640	4.0						STEA7640T	
	EZ Post 7650	5.0						STEA7650T	
	EZ Post 7610	1.0	Ø7.0	5.5	Non-Hex		STEA7610NT		
	EZ Post 7620	2.0					STEA7620NT		
	EZ Post 7630	3.0					STEA7630NT		
	EZ Post 7640	4.0					STEA7640NT		
	EZ Post 7650	5.0					STEA7650NT		

# Multi Unit Abutment



\* **EZ Post Cylinder** In the case of MUA S Type, the EZ Post Cylinder works like a scan abutment.

# Scan abutment

Scan Abutment **MUA-N Type**



Product info			Library info (Choose below CAD library option)		
Diameter	Lenght	Ref.C	exocad library name	3shape library name	
Ø4.0	13	AMUASR4013T	MUA Custom Abutment (13mm)	MUA_ZrGEN_Custom abutment	

Scan Abutment **MUA-S Type**



Product info			Library info (Choose below CAD library option)		
Diameter	Lenght	Ref.C	exocad library name	3shape library name	
Ø5.0	9	RCA900T	MUA S_Type_EZpost-Cylinder	-	


## ZrGEN

Ti-Base



Library info		Product info				
ZrGEN Option In exocad library		Type	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Ref.C
System	Kit					
 MUA Custom Abutment (13mm) MUA-N Type	ZrGEN 5515	N-Type	5.5	0.8	5	AMUAPR5515N.MTN
	ZrGEN 5516				6	AMUAPR5516N.MTN
	ZrGEN 5518				8	AMUAPR5518N.MTN

## Bar



Library info		Product info				
Bar Option In exocad library		Type	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Ref.C
System	Kit					
 MUA N Type Bar MUA-N Type	Bar	N-Type	-	-	-	-

## EZ Post Cylinder



Library info		Product info			
ZrGEN Option In exocad library		Type	Diameter (Ø)	Post Height (mm)	Ref.C
System	Kit				
 MUA S Type_EZpost Cylinder MUA-S Type	EZ Post Cylinder	S-Type	5.0	5	RCA900T





**ZrGEN**

Ti-Base




MUA ZrGEN  
Custom Abutment

**MUA-N Type**

Library info		Product info				
ZrGEN Option In 3shape library		Type	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Ref.C
System	Kit for 13mm					
	ZrGEN 5515_S13	N Type	5.5	0.8	5	AMUAPR5515N.MTN
	ZrGEN 5516_S13				6	AMUAPR5516N.MTN
	ZrGEN 5518_S13				8	AMUAPR5518N.MTN

# Octa



# Scan abutment

## Scan Abutment [Octa]



Product info			Library info (Choose below CAD library option)	
Diameter	Lenght	Ref.C	exocad library name	3shape library name
Ø4.0	11	AOCEC4011T	Octa_Level_Custom Abutment (11mm)	Octa Level_ZrGEN_Custom abutment Octa Level_TiGEN_Custom abutment



**ZrGEN**

Ti-Base



Octa Level  
Custom Abutment (11mm)

Library info		Product info					
ZrGEN Option In exocad library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit						
	ZrGEN 5015	Small	Ø5.0	0.8	5.0	Octa	AOCEPS5015.MTN
	ZrGEN 5515	Regular	Ø5.5				ANOEPR5515.MTN
	ZrGEN 6515	Wide	Ø6.5				AOCEPS6515.MTN
	ZrGEN 5015	Small	Ø5.0	0.8	5.0	Non-Octa	ANOEPS5015.MTN
	ZrGEN 5515	Regular	Ø5.5				ANOEPR5515.MTN
	ZrGEN 6515	Wide	Ø6.5				ANOEPW6515.MTN



**TiGEN**  
Pre-milled




Library info		Product info				
Pre-milled Option In exocad library		Size	Diameter (∅)	Length (mm)	Type	Ref.C
System	Kit					
<p>Octa Level Custom Abutment (11 mm)</p>	TiGEN 1020	Small	10	20	Hex	OCTS1020.MTN
	TiGEN 1220		12	20		OCTS1220.MTN
	TiGEN 1020	Regular	10	20		OCTR1020.MTN
	TiGEN 1220		12	20		OCTR1220.MTN
	TiGEN 1020	Wide	10	20		OCTW1020.MTN
	TiGEN 1220		12	20		OCTW1220.MTN
	TiGEN 1020	Small	10	20	Non-Hex	NOTS1020.MTN
	TiGEN 1220		12	20		NOTS1220.MTN
	TiGEN 1020	Regular	10	20		NOTR1020.MTN
	TiGEN 1220		12	20		NOTR1220.MTN
	TiGEN 1020	Wide	10	20		NOTW1020.MTN
	TiGEN 1220		12	20		NOTW1220.MTN



**ZrGEN**

Ti-Base




Library info		Product info					
ZrGEN Option In 3shape library		Size	Diameter (Ø)	Cuff Height (mm)	Post Height (mm)	Type	Ref.C
System	Kit for 11mm						
 <p>Octa Level ZrGEN Custom Abutment (11 mm)</p>	ZrGEN 5015_S11	Small	Ø5.0	0.8	5.0	Octa	AOCEPS5015.MTN
	ZrGEN 5515_S11	Regular	Ø5.5				ANOEPS5515.MTN
	ZrGEN 6515_S11	Wide	Ø6.5				ANOEPS6515.MTN
	ZrGEN 5015_S11	Small	Ø5.0	0.8	5.0	Non-Octa	ANOEPS5015.MTN
	ZrGEN 5515_S11	Regular	Ø5.5				ANOEPR5515.MTN
	ZrGEN 6515_S11	Wide	Ø6.5				ANOEPW6515.MTN



**TiGEN**  
Pre-milled



Library info		Product info				
Pre-milled Option In 3shape library		Size	Diameter (Ø)	Length (mm)	Type	Ref.C
System	Kit for 11mm					
 <p>Octa Level TiGEN Custom Abutment (11 mm)</p>	TiGEN 1020	Small	10	20	Hex	OCTS1020.MTN
	TiGEN 1220		12	20		OCTS1220.MTN
	TiGEN 1020	Regular	10	20		OCTRI020.MTN
	TiGEN 1220		12	20		OCTRI220.MTN
	TiGEN 1020	Wide	10	20		OCTW1020.MTN
	TiGEN 1220		12	20		OCTW1220.MTN
	TiGEN 1020	Small	10	20	Non-Hex	NOTS1020.MTN
	TiGEN 1220		12	20		NOTS1220.MTN
	TiGEN 1020	Regular	10	20		NOTRI020.MTN
	TiGEN 1220		12	20		NOTRI220.MTN
	TiGEN 1020	Wide	10	20		NOTW1020.MTN
	TiGEN 1220		12	20		NOTW1220.MTN





**Contact us**

[r2b@r2gate.com](mailto:r2b@r2gate.com)  
[www.r2gate.com](http://www.r2gate.com)

Made by Simon  
[r2.sales3@imegagen.com](mailto:r2.sales3@imegagen.com)